


MICROGEN BIOPRODUCTS LTD



MICROGEN[®] STREP-ID

AN IDENTIFICATION SYSTEM FOR STREPTOCOCCI AND
ENTEROCOCCI FROM CLINICAL, VETERINARY, FOOD AND
ENVIRONMENTAL SOURCES.

PRODUCT CODE: MID-62 (20 TEST KIT)

- ◆ Substrates selected specifically for Streptococcus, Enterococcus and related species
- ◆ Simple and easy-to-use
- ◆ Results in 24 hours
- ◆ Results interpreted using new Microgen ID system software (registered users are entitled to free database updates)
- ◆ Excellent differentiation of *Enterococci* spp. especially *E. casseliflavus*, *E. mundtii* and *E. gallinarum*

Protecting Food and Health

MICROGEN BIOPRODUCTS LTD

1 Admiralty Way
Camberley
Surrey
GU15 3DT
United Kingdom

Phone: +44 1276 600081
Fax: +44 1276 600151
E-mail: sales@microgenbioproducts.com
www.microgenbioproducts.com



Microgen® Strep-ID is a biochemical test system which utilises a 12 well (12 test) microwell test strip and 3 off-strip tests; hippurate hydrolysis (provided), alpha-haemolysis and beta-haemolysis for the identification of Streptococcal and Enterococcal species important to clinical, veterinary, food and environmental laboratories.

Performance Characteristics:

Microgen® Strep-ID has been evaluated in comparison with a well-established commercially available product for identification of cultured bacterial isolates. A total of 69 cultures, comprising cultures from recognised culture collections and clinical isolates were examined.

	Number Tested	MID-62	Competitor
<i>E. avium</i>	1	1	1
<i>E. durans</i>	2	2	1
<i>E. faecalis</i>	10	10	10
<i>E. faecium</i>	5	5	5
<i>E. gallinarum</i>	1	1	0
<i>E. hirae</i>	1	1	0
<i>G. haemolysans</i>	1	1	1
<i>S. acidominimus</i>	1	1	0
<i>S. agalactiae</i>	5	5	5
<i>S. anginosus</i>	3	3	2
<i>S. bovis</i>	1	1	1
<i>S. constellatus</i>	3	2	3
<i>S. dysgalactiae ssp. equi</i>	1	1	1
<i>S. equi subsp. equi</i>	2	2	2
<i>S. equi subsp. zooepidemicus</i>	1	1	1
<i>S. gordonii</i>	1	1	0
<i>S. intermedius</i>	1	1	0
<i>S. mitis</i>	9	9	5
<i>S. mutans</i>	4	3	4
<i>S. parva</i>	1	1	0
<i>S. pneumoniae</i>	5	4	2
<i>S. pyogenes</i>	3	3	3
<i>S. salivarius</i>	4	4	3
<i>S. sanguinis</i>	2	0	2
<i>S. uberis</i>	1	1	1
<i>S. vestibularius</i>	1	1	0
TOTAL	70	65	53

The Microgen® Strep-ID correctly identified 70 (93%) and the competitor product correctly identified 52 (74%) of the isolates tested. The superior specification of Microgen® Strep-ID is complimented by an excellent differentiation of Enterococci spp. and a reduced test system minimises handling and manipulation whilst reducing procedural error and time.

As with all Microgen Bioproducts biochemical identification systems Microgen Strep-ID is supported by the simplest-to-use and most feature-packed computer identification system (product code: MID60) available, delivering the most comprehensive and up-to-date databases which registered users can update free-of-charge from www.microgenbioproducts.com

Microgen® Strep-ID:

Code: MID-62
Description: Microgen® Strep-ID
Size: 20 Test
Storage: 2-8°C

Additional Reagents Available:

MID61A Nitrate A reagent
MID61B Nitrate B reagent
MID61C VP I reagent
MID61D VP II reagent
MID61F Indole Kovacs reagent
MID61G Oxidase Strips
MID61H Mineral oil
MID61K PYR reagent
MID61L Ninhydrin reagent

Other Microgen® Identification Systems Available:

MID-64 Microgen® GN A –ID
MID-65 Microgen® GN B - ID
MID-66 Microgen® Bacillus - ID
MID-67 Microgen® Listeria - ID
MID-69 Microgen® Staph – ID