



**PROCESS FILTRATION**  
**PRODUCT GUIDE**

# Company Profile

Advanced Microdevices (**mdi**) is a leader in innovative membrane technologies. Starting from a single person R&D operation in 1976, **mdi** has developed into a dedicated team of 300 plus, manufacturing more than 15000 products.

The company's core competence is its ability to develop new membrane technologies and innovate existing ones to deliver advantages to the customer for high end purification and separation applications in a wide range of industries such as pharmaceuticals, biopharmaceuticals, biotechnology, food and beverage, hospitals, and immunodiagnostics.



**Existing Facility (100,000 sq. ft.)**



**New GMP Facility (100,000 sq. ft.)**



**mdi** is a vertically integrated company with inter disciplinary expertise that helps deliver an ever increasing portfolio of innovative products with unique performance advantages.

**mdi** products are used in highly regulated industries for critical applications such as sterile filtration of pharmaceutical and biopharmaceutical process fluids, sterility testing of aseptic drug products, analytical sample preparation, new drug development, and making reliable immunoassays for testing of diseases at patient bedside.

**mdi's** modern GMP facilities with large ISO 7 Clean Areas more than meets the required standards. The products are manufactured in ISO 9001:2008 certified facility with the help of trained manpower meeting or exceeding industry standards. Many **mdi** products are recognized as the best available in the world.

# Table of Contents

|  |                 |  |  |       |
|--|-----------------|--|--|-------|
| <b>General</b>   | <b>Page No.</b> |  | <b>Filters for Clarification and pre-filtration</b>          | 39    |
| ➤ MembraneTechnologies   | 2               |  | ➤ Filter Selection Chart                                     | 40    |
| ➤ Quality Assurance  | 3               |  | ➤ ClariPro GK- 1", 2", 5", 8"                                | 41    |
| ➤ Validation Services  | 4               |  | ➤ ClariPro GK- 10", 20", 30"                                 | 42    |
| ➤ Customer Support   | 5               |  | ➤ ClariSure GK   | 43    |
| ➤ Filter Selection and Sizing  | 6               |  | ➤ ClariCap GS- 1", 2", 5", 8"                                | 44    |
| ➤ Filter Sizing: Linear Upscaling from R&D to Production Process                   | 7               |  | ➤ ClariCap GS- 10", 20", 30"                                 | 45    |
| ➤ Easy Connect   | 8               |  | ➤ ClariSure GS   | 46    |
| <b>Filters for Sterilization of Liquids: PES Membrane Filters</b>                  | 9               |  | ➤ ClariSure GP   | 47    |
| ➤ Filter Selection Chart   | 10              |  | ➤ ClariCap PP- 1", 2", 5", 8"                                | 48    |
| ➤ AseptiCap KL/KS- 25mm, 50mm  | 11              |  | ➤ ClariCap PP- 10", 20", 30"                                 | 49    |
| ➤ AseptiCap KL/KS- 1", 2", 5", 8"  | 13              |  | ➤ ClariSure PA   | 50    |
| ➤ AseptiCap KS- 10", 20", 30"  | 15              |  | ➤ Microglassfiber Disc Filters                               | 51    |
| ➤ AseptiCap KO - 1", 2", 5", 8"  | 17              |  | <b>Filters for Polishing and Clarification</b>               | 52    |
| ➤ AseptiSure HS  | 18              |  | ➤ Filter Selection Chart                                     | 53    |
| ➤ AseptiSure KS  | 19              |  | ➤ ClariSure PP   | 54    |
| ➤ AseptiSure KR  | 20              |  | ➤ ClariSure PL   | 55    |
| <b>Filters for Sterilization of Liquids Nylon-66 Membrane Filters</b>              | 21              |  | ➤ ClariSure DP   | 56    |
| ➤ Filter Selection Chart   | 22              |  | <b>Chemical Compatibility</b>                                | 57-58 |
| ➤ AseptiCap NL/NS- 1", 2", 5", 8"  | 23              |  | <b>End Connection Availability Chart for Capsule Filters</b> | 59    |
| ➤ AseptiCap NS- 10", 20", 30"  | 24              |  | <b>Dimensions Chart for Capsule Filters</b>                  | 59    |
| ➤ AseptiSure NS  | 25              |  | <b>Ordering Information</b>                                  | 60    |
| ➤ Nylon-66 Membrane Disc Filters   | 26              |  |  |       |
| <b>Filters for Sterilization of Air/Gas</b>  | 27              |  |  |       |
| ➤ Filter Selection Chart   | 28              |  |  |       |
| <b>Single Use Gamma Sterilizable PVDF Membrane Filter Devices</b>                  |                 |  |  |       |
| ➤ AseptiVent VF - 25mm, 37mm, 50mm   | 29              |  |  |       |
| ➤ AseptiVent VF - 1", 2", 5", 8"   | 31              |  |  |       |
| <b>Multiple Use Autoclavable / Steam Sterilizable PTFE Membrane Filter Devices</b> |                 |  |  |       |
| ➤ AseptiVent TF- 25mm, 37mm, 50mm  | 33              |  |  |       |
| ➤ AseptiVent TF- 1", 2", 5", 8"  | 35              |  |  |       |
| ➤ AseptiVent TF- 10"   | 36              |  |  |       |
| ➤ AseptiSure TH  | 37              |  |  |       |
| ➤ AseptiSure TF  | 38              |  |  |       |

mdi offers customized solutions to enhance process efficiency, productivity, product quality, and consistency in wide ranging applications in biopharmaceuticals, pharmaceuticals, biotechnology, microelectronics, and immunodiagnostics.

## Research and Development

A unique multidimensional research and development facility at mdi, involving teams specializing in critical areas of membrane technology, biotechnology, electronics, chemistry, and mechanical engineering, continuously strives with an integrated approach to develop advanced, innovative, and customized products.

These filters have innovative design inputs to deliver unique performance advantages over competing products in terms of higher retention efficiencies, flow rates, throughputs, and lower filtration losses.

This unparalleled capability to develop custom products and solutions is a continuous source of pride and drives the mdi team to push the boundaries of technology and maximize value for the user.



mdi Microfiltration products are well designed with built-in quality assurance.

The careful selection of raw materials, validated production processes and Quality Management System certified by ISO 9001:2008 ensures consistently high quality products. mdi products meet 21 CFR, ASTM, compendia requirements and meet global regulatory expectations.

## Facilities

mdi filters are produced and packaged in facilities meeting GMP requirements such as Clean Rooms with Class 10,000 and Class 100 areas for critical processes.

## Deep Characterization and Certification

Apart from retention efficiency and other functional parameters such as flowrates, temperature/hydraulic stress testing etc., mdi filters are deeply characterized for critical areas of concern such as biosafety, bioburden levels, endotoxin levels, and extractables.

## Quality Control

The filters go through stringent in process and final product testing and quality is ensured by in place QMS.

## Traceability

Each sterilizing grade cartridge and capsule filter has unique identification number and is accompanied with individual certificate of quality.



As per regulatory requirement, the pharmaceutical industry has to provide a high level of assurance that the sterile drug product manufactured through aseptic processing offers the identity, strength, quality, and purity it purports to have or is represented to possess (Ref. USFDA 21CFR 211.100(a)). Consequently it has become increasingly critical to establish/quantify the impact on the drug due to its interface with various process components under different process conditions.

Sterilizing grade filters are of critical importance in aseptic manufacturing, and it is the drug manufacturer's responsibility to show that the selected filter is able to sterilize the product under the process conditions and it does not affect the purity, quality and strength of the drug product.

## Validation Services

**mdi ascertain** Filter validation services are designed to meet customer specific needs and help achieve regulatory compliance. These include the following:

- **Studies establishing filter integrity test values specific to drug product**
- **Filter fluid interaction studies**
  - Physico-Chemical compatibility studies
  - Extractable/Leachable studies
  - Adsorption studies
- **Microbial retention studies**

All of these studies are executed as per pre-approved test methodologies establishing the test conditions and acceptance criteria.

**mdi** also offers post validation support for regulatory audits.

## Validation Guides

**mdi** filters are validated as per global regulatory requirements. These filters, in compliance with the Regulation Title 21 Code of Federal Regulations (CFR) Part 314.420, have been registered at the U.S. Food and Drug Administration through Drug Master File # DMF 15554

Detailed documentation on validation of **mdi** filters for sterilization of fluids (air/gases and liquids) in form of Validation Guides is also available.



mdi technology executives assist in problem solving and process upgradation through experience sharing and developing customized products and systems. Some of these customer support activities are:

## » Customized Filtration Solutions

mdi offers customized solutions for complex filtration problems. mdi technology executives will help you in finding solutions to filter difficult to filter fluids & minimizing filtration losses.

## » Filtration System Design

Designing an efficient filtration system is an integral part of process optimization for minimizing filtration costs, increasing yields, and reducing process time. mdi offers technical support for selection of filter materials, performing throughput studies to optimize filter train and filter sizing.

## » Installation and Operation Qualification

mdi provides well documented installation, operation and performance qualification guidelines for all equipment and systems it offers.

## » Regulatory Assistance

mdi provides complete regulatory assistance to its customers. mdi products and validation services meet global regulatory expectations.

## » Technical Seminars

Technical seminars at customer's location are organized to help the users understand the performance characterization of product in use, differentiate between various options available in the market and select the best solutions to suit their requirements. These interactions help create optimized systems and also upgrade current processes in terms of performance and cost effectiveness.



# Filter Selection and Sizing

Highly regulated process industries such as pharmaceuticals and biopharmaceuticals are working with a very wide variety of process streams/fluids under different process conditions. They are continuously faced with the challenge of achieving their process objectives efficiently and cost effectively.

Microfiltration is a key process step to achieve critical process objectives that range from sterilization of process fluids to bioburden reduction, polishing and clarification.

These process streams can range from easy to filter water for injection to difficult to filter colloidal solutions, emulsions, liposomal drug delivery systems or large molecule high value therapeutic proteins and vaccine concentrates. Such a wide spectrum of process streams coupled with different process objectives is quite a challenge for the process owner.

Selection of the right filters and their sizing to meet various process needs is thus critical to successfully achieve the desired objectives.

This would require the following important questions to be answered with respect to the process:

1. What is the objective of filtration?
  - Sterilization
  - Bioburden reduction
  - Particle removal
  - Clarification
2. What is the fluid to be filtered?
  - Liquid or gases
  - What are the fluid properties such as pH, viscosity, temperature and surface tension?
3. What will be the process conditions?
  - Will the filtration system be inline steam sterilized or autoclaved?
  - Will the system be sanitized with chemicals or hot water?
  - What are the sterilization/sanitization conditions?
  - Are the filters going to be used once or multiple times?
  - How many times the system will be sterilized/sanitized?
  - What is the maximum operating temperature?
  - What is the allowable maximum differential pressure?

4. Is it going to be a continuous process or batch filtration?
5. What will be the batch volume for full scale process filtration?
6. What is the maximum allowable filtration time or the desired minimum flow rate?

Once the filter type with respect to MOC, pore size etc. has been established the next step is to establish the filter size.

Filter sizing, although to some extent is dependent on factors such as minimum desired flow rate, fluid viscosity as well as temperature, a critical parameter in the contamination profile of the fluid to be filtered. The nature and quantum of contaminants defines their interaction with the filter, which in turn defines the throughput one can achieve from a given filter area for the fluid in question. An understanding of this behavior will help define not only the right sized filtration system but also the right combination of pre-filters to achieve desired/optimum throughputs.

**mdi** offers filter sizing services to product development labs and process owners in full scale manufacturing. These involve small scale throughout studies to establish suitable and cost effective filtration system. Different lab scale filter and pre-filter combinations are used to maximize throughputs. The selected combination is, based on desired batch volumes or throughputs, linearly extrapolated to establish filter size.

**For more information please contact our local technology executive or write to [support@mdimembrane.com](mailto:support@mdimembrane.com)**



# Filter Sizing: Linear Upscaling from R&D to Production Process

Researchers in NDDS and formulation development are concerned about filter fluid interaction impacting the stability, purity, strength etc. of the drug product, and they take a keen interest in filter selection at the formulation development stage itself.

Although preliminary compatibility data support initial filter selection, for stability studies detailed filter validations are required to provide enough documented evidence to justify specific filter use.

A critical requirement that needs to be addressed at this stage is of scalability from R&D to pilot scale to full scale production processes. Any change in filter MOC for full scale processes will require additional validation.

**mdi** offers a wide range of filters to provide linear scale up from lab scale to production process. While scaling up the process, the appropriate filter size can be selected by increasing the effective filtration area of filter proportionate to the process fluid volumes.

All Materials of construction of core, sleeve, end caps, support layers and housing as well as manufacturing process is identical for all filter devices starting from 5 cm<sup>2</sup> to 18000 cm<sup>2</sup> hence process scaling can be facilitated without triggering additional validation studies for given process conditions. **mdi** provides complete documentation for each of the filters thereby reducing the additional validation cost and time.



25 mm, 5 cm<sup>2</sup>



50 mm, 20 cm<sup>2</sup>



1", 250 cm<sup>2</sup>



2", 500 cm<sup>2</sup>



5", 1000 cm<sup>2</sup>



8", 2000 cm<sup>2</sup>



10", 6000 cm<sup>2</sup>

# Easy Connect

## Wide Range of End Connections

Pharmaceutical and Biopharmaceutical processes involve transfer of high value fluids through multiple process steps. Making high quality, reliable, flexible and functionally convenient connectivity with filters is of utmost value.

**mdi** Capsule filters offer a wide range of reliable end connections for functional convenience and customized connectivity.

## Validated for Performance

These end connections are manufactured with tight dimension tolerance and are validated for strength and connection integrity under extreme use and sterilization conditions.

## Customized Connectivity

**mdi** Capsule filters can also be customized to offer different inlet-outlet combinations to meet the unique connectivity needs in process assemblies where, for example, stainless steel components with sanitary flange connections are sometimes required to be connected to single use disposable systems through quick-connectors or hose barb connections.



**3/4" Sanitary Flange**



**1 1/2" Sanitary Flange**



**1/2" HB**



**Single Stepped Hose Barb**



**1/4" SHB**



**Quick Connector**



**1 1/2" Sanitary Flange to 1/2" Barb Hose**



**1 1/2" Sanitary Flange to 3/4" Sanitary Flange**



**AseptiCap with HighSecurity 1/2" hose barb connection**

# Filters for Sterilization of Liquids: PES Membrane Filters

**mdi** produces a wide range of Gamma sterilizable and steam sterilizable PES membrane capsule and cartridge filters to meet filtration requirements of biopharmaceutical and pharmaceutical processing.

These filters meet key process requirements such as high retention efficiency, very high protein recoveries, extremely low extractables, high throughputs, wide chemical compatibility etc.

**mdi** PES filter devices are available as:

| Filter Type   | Single Layer                               | Multiple Layer       |
|---|--|----------------------|
| Gamma Sterilizable Capsule Filters                              | <i>AseptiCap KL</i>                        | <i>AseptiCap KS</i>  |
| Autoclavable Capsule Filters                                    | <i>AseptiCap KL</i><br><i>AseptiCap KO</i> | <i>AseptiCap KS</i>  |
| Steam Sterilizable Cartridge Filters                            | <i>AseptiSure KR</i>                       | <i>AseptiSure KS</i> |
| High Temperature Resistant Steam Sterilizable Cartridge Filters | –  | <i>AseptiSure HS</i> |

## Quality Assurance

These filter devices are manufactured in class 10,000 clean rooms under ISO 9001 : 2008 certified quality management systems and are validated to meet compendia and regulatory requirements.

## Applications

### Sterile Filtration of:

- Cell culture media
- Cell culture media containing serum
- Media additives
- Final product concentrates
- Bioburden reduction from cell harvest supernatants
- Filtration of equilibrating, washing and elution buffers for chromatography steps
- Sterilization of biopharmaceuticals such as vaccines and therapeutic proteins
- Sterilization of buffers
- Sterilization of adjuvants
- Manufacture of high value SVP like Oncology drugs, hormones and therapeutic proteins
- Manufacture of Large Volume Parenterals
- Manufacture of Water for Injection
- Manufacture of Small Volume Parenterals

| Assurance                          |  |
|------------------------------------|--|
| <b>Toxicity</b>                    | Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics  |
| <b>Cytotoxicity</b>                | Passes Biological Reactivity Tests, In Vitro, USP <87> for cytotoxicity  |
| <b>Bioburden</b>                   | Bioburden level is < 1000 cfu/filter device as per ANSI/AAMI/ISO 11737-1 : 1995  |
| <b>Bacterial Endotoxin</b>         | Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test  |
| <b>Non Fiber Releasing</b>         | Passes test as per USP and comply with USFDA 21 CFR Part 211.72 for fiber release  |
| <b>Extractables with WFI</b>       | Passes test as per USP   |
| <b>Oxidizable Substances</b>       | Within limits as specified in USP  |
| <b>Particle Shedding</b>           | Passes USP test for particulates in injectables  |
| <b>TOC/Conductivity at 25 °C</b>   | Meets the WFI requirements of USP <643> for Total Organic Carbon and USP <645> for Water Conductivity after a specified volume of purified water flush |
| <b>Indirect Food Additive</b>      | All Polypropylene components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177.1520   |
| <b>Good Manufacturing Practice</b> | These products are manufactured in a facility which adheres to Good Manufacturing Practices  |

# Filter Selection Chart

| Application Area   | Key Application Requirements   | Gamma Sterilizable Capsule Filters                                  | Steam Sterilizable  |  |   |
|--|--|---|---|--|---|
|  |  |   | Capsule Filters   | Cartridge Filters  |   |
| <b>Biopharmaceuticals</b>                                      |  |   |   |  |   |
| Media preparation  | - Mycoplasma removal (in case of Mammalian Cell Culture)                       | <b>AseptiCap KS</b> - 0.1 µm PES Membrane Capsule Filter            | <b>AseptiCap KS</b> 0.1 µm PES Membrane Capsule Filter            | <b>AseptiSure KS</b> 0.1 µm PES Membrane Cartridge Filter            | <b>AseptiSure HS</b> 0.1 µm High Temperature PES Membrane Cartridge Filter            |
|  | - Microbial retention (in case of Microbial Fermentation)                      | <b>AseptiCap KS</b> - 0.2 µm PES Membrane Capsule Filter            | <b>AseptiCap KS</b> 0.2 µm PES Membrane Capsule Filter            | <b>AseptiSure KS</b> 0.2 µm PES Membrane Cartridge Filter            | <b>AseptiSure HS</b> 0.2 µm High Temperature PES Membrane Cartridge Filter            |
| Sterile filtration of growth regulators                        | - Absolute retention<br>- Low protein binding                                  | <b>AseptiCap KS</b> - 0.2 µm PES Membrane Capsule Filter            | <b>AseptiCap KS</b> 0.2 µm PES Membrane Capsule Filter            | <b>AseptiSure KS</b> 0.2 µm PES Membrane Cartridge Filter            | <b>AseptiSure HS</b> 0.2 µm High Temperature PES Membrane Cartridge Filter            |
| Sterile filtration of alkaline/acidic solutions for pH control | - Absolute retention<br>- Compatible with 1-14 pH                              | -   | <b>AseptiCap KO</b> 0.2 µm PES Membrane Capsule Filter            | <b>AseptiSure KR</b> 0.2 µm PES Membrane Cartridge Filter            | -   |
| Cell Harvesting  | - Bioburden reduction  | <b>AseptiCap KS</b> - 0.2 µm or 0.45 µm PES Membrane Capsule Filter | <b>AseptiCap KS</b> 0.2 µm or 0.45 µm PES Membrane Capsule Filter | <b>AseptiSure KS</b> 0.2 µm or 0.45 µm PES Membrane Cartridge Filter | <b>AseptiSure HS</b> 0.2 µm or 0.45 µm High Temperature PES Membrane Cartridge Filter |
| Buffer filtration  | - Bioburden reduction<br>- Absolute retention                                  | <b>AseptiCap KS</b> - 0.2 µm or 0.45 µm PES Membrane Capsule Filter | <b>AseptiCap KS</b> 0.2 µm or 0.45 µm PES Membrane Capsule Filter | <b>AseptiSure KS</b> 0.2 µm or 0.45 µm PES Membrane Cartridge Filter | <b>AseptiSure HS</b> 0.2 µm or 0.45 µm High Temperature PES Membrane Cartridge Filter |
| Sterile filtration of vaccines and therapeutic proteins        | - Absolute retention<br>- Low protein binding<br>- Low holdup volume           | <b>AseptiCap KS</b> - 0.2 µm PES Membrane Capsule Filter            | <b>AseptiCap KS</b> 0.2 µm PES Membrane Capsule Filter            | -  | -   |
| <b>Pharmaceuticals</b>   |  |   |   |  |   |
| Large Volume Parenterals                                       | - Absolute retention<br>- High throughputs                                     | -   | -   | <b>AseptiSure KS</b> 0.2 µm PES Membrane Cartridge Filter            | <b>AseptiSure HS</b> 0.2 µm High Temperature PES Membrane Cartridge Filter            |
| Small Volume Parenterals                                       | - Absolute retention<br>- Low Protein Binding<br>- Wide Chemical Compatibility | <b>AseptiCap KS</b> - 0.2 µm PES Membrane Capsule Filter            | <b>AseptiCap KS</b> 0.2 µm PES Membrane Capsule Filter            | <b>AseptiSure KS</b> 0.2 µm PES Membrane Cartridge Filter            | <b>AseptiSure HS</b> 0.2 µm High Temperature PES Membrane Cartridge Filter            |
| WFI  | - Absolute retention   | <b>AseptiCap KS</b> - 0.2 µm PES Membrane Capsule Filter            | <b>AseptiCap KS</b> 0.2 µm PES Membrane Capsule Filter            | <b>AseptiSure KS</b> 0.2 µm PES Membrane Cartridge Filter            | <b>AseptiSure HS</b> 0.2 µm High Temperature PES Membrane Cartridge Filter            |

## AseptiCap KL/KS- 25 mm and 50 mm

**Microbially Validated  
as per ASTM F 838-05**

**Complies with  
USFDA 21 CFR 211.72**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

Specially designed filters for process development and formulation development labs with identical materials of construction for easy scale up to larger capsule filters and cartridge filters.

50mm is a specially vented device for use with peristaltic pump to ensure easy removal of entrapped air in the upstream.

**Radiation Sterilizable:** AseptiCap KL/KS-

**Autoclavable:** AseptiCap KL/KS




### Specifications


| Construction                                     |   |   |  |
|--|---|---|--|
| Final Filter Pore Size                           | 0.1 µm  | 0.2 µm  | 0.45 µm  |
| Prefilter Pore Size<br>(in case of AseptiCap KS) | 0.45 µm   | 0.8 µm, 0.65 µm, 0.45 µm  | 0.8 µm, 0.65 µm  |
| Membrane   | Hydrophilic PES   |   |  |
| Support Layers                                   | Polyester   |   |  |
| Body and Core                                    | Polypropylene   |   |  |
| Integrity Testing/Retention                      |   |   |  |
| Bubble Point                                     | ≥ 31 psi (2.18 Kg/cm <sup>2</sup> )<br>with 50% IPA/Water Solution              | ≥ 50 psi (3.51 Kg/cm <sup>2</sup> )<br>with Water   | ≥ 30 psi (2.11 Kg/cm <sup>2</sup> )<br>with Water                  |
| Microbial Retention                              | LRV >7 for Acholeplasma laidlawii<br>(ATCC 23206) per cm <sup>2</sup>           | LRV >7 for Brevundimonas diminuta<br>(ATCC 19146) per cm <sup>2</sup>                     | LRV >7 for Serratia marcescens<br>(ATCC 14756) per cm <sup>2</sup> |
| Size   |   |   |  |
| Size   | 25 mm   | 50 mm   |  |
| Effective Filtration Area (Nominal)              | 5 cm <sup>2</sup>   | 20 cm <sup>2</sup>  |  |
| Operational                                      |   |   |  |
| Max. Operating Temperature                       | 55 °C   | 60 °C   |  |
| Max. Differential Pressure                       | 75 Psi (5 Kg/cm <sup>2</sup> )@ 25 °C   | 42 Psi (3 Kg/cm <sup>2</sup> ) @ 30 °C  |  |
| Hold-up Volume(with air purge)                   | <50 µL  | <300 µL   |  |
| Burst Pressure                                   | > 14 Kg/cm <sup>2</sup>   | > 8 Kg/cm <sup>2</sup>  |  |
| Sterilization                                    | By Irradiation  | AseptiCap KL/KS - Gamma Irradiatable up to 50 kGy   |  |
|  | By Gas  | AseptiCap KL/KS: Sterilizable by Ethylene Oxide   |  |
|  | By Autoclave  | AseptiCap KL/KS - : Autoclavable at 121°C for 30 minutes, 1 cycle after gamma irradiation |  |
|  |   | AseptiCap KL/KS: Autoclavable at 121°C for 30 minutes, 25 cycles                          |  |
|  |   | <b>These cannot be inline steam sterilized</b>  |  |
| Shelf Life                                       | 2 years after gamma sterilization<br>3 years after Ethylene Oxide sterilization |   |  |
| pH Compatibility                                 | Compatible with pH range of 1-10  |   |  |

## Water Flow Rates

AseptiCap KL - 25 mm



AseptiCap KL - 50 mm



## Ordering Information

AseptiCap KL/KS and AseptiCap KL/KS - , 25 mm

| Type                            |      | Size  |      | Pore Size |      | Inlet/Outlet     |      | Radiation Sterilizable |      | X | Sterility     |      | Pack Size |      |
|---------------------------------|------|-------|------|-----------|------|------------------|------|------------------------|------|---|---------------|------|-----------|------|
|                                 | Code | Dia   | Code |           | Code |                  | Code |                        | Code |   |               | Code |           | Code |
| AseptiCap KL                    | IKLX | 25 mm | 06   | 0.1 µm*   | 36   | Female Luer Lock | M    | Yes                    | R    | X | Non Sterile   | 1    | 100       | 04   |
| AseptiCap KS (0.45 µm Upstream) | IKSX |       |      | 0.2 µm    | 01   | Male Luer Slip   | N    | No                     | X    |   | EO Sterile    | 2    |           |      |
| AseptiCap KS (0.65 µm Upstream) | IKS3 |       |      | 0.45 µm   | 02   |                  |      |                        |      |   | Gamma Sterile | 3    |           |      |
| AseptiCap KS (0.8 µm Upstream)  | IKS5 |       |      |           |      |                  |      |                        |      |   |               |      |           |      |
| <b>Example</b>                  |      |       |      |           |      |                  |      |                        |      |   |               |      |           |      |
| IKLX                            |      | 06    |      | 01        |      | MN               |      | R                      |      | X | 1             |      | 04        |      |

AseptiCap KL/KS and AseptiCap KL/KS - , 50 mm

| Type                            |      | Size  |      | Pore Size |      | Inlet/Outlet       |      | Radiation Sterilizable |      | Bell | Sterility |               | Pack Size |    |    |
|---------------------------------|------|-------|------|-----------|------|--------------------|------|------------------------|------|------|-----------|---------------|-----------|----|----|
|                                 | Code | Dia   | Code |           | Code |                    | Code |                        | Code |      | Code      |               | Code      |    |    |
| AseptiCap KL                    | VKLX | 50 mm | 10   | 0.1 µm*   | 36   | ¼" SHB             | B    | Yes                    | R    | Yes  | B         | Non Sterile   | 1         | 12 | 08 |
| AseptiCap KS (0.45 µm Upstream) | VKSX |       |      | 0.2 µm    | 01   | ¾" Sanitary Flange | S    | No                     | X    | No   | X         | EO Sterile    | 2         |    |    |
| AseptiCap KS (0.65 µm Upstream) | VKS3 |       |      | 0.45 µm   | 02   |                    |      |                        |      |      |           | Gamma Sterile | 3         |    |    |
| AseptiCap KS (0.8 µm Upstream)  | VKS5 |       |      |           |      |                    |      |                        |      |      |           |               |           |    |    |
| <b>Example</b>                  |      |       |      |           |      |                    |      |                        |      |      |           |               |           |    |    |
| VKSX                            |      | 10    |      | 36        |      | SS                 |      | X                      |      | X    | 1         |               | 08        |    |    |

\* Note: 0.1 µm is available in AseptiCap KS and AseptiCap KS - only with 0.45 µm upstream

For End Connection availability and dimensions with different sizes refer Page no. 59.

Note: Gamma Sterile filters can not be gamma irradiated again

Example for Non Sterile: VKSX1036SSRX104

Example for gamma Sterile: VKSX1036SSXX304

# AseptiCap KL/KS- 1", 2", 5", 8"

**Microbially Validated  
as per ASTM F 838-05**

**Complies with  
USFDA 21 CFR 211.72**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

Polyethersulfone membrane capsule filters are self contained, ready to use, disposable filtration devices that contain a mini cartridge filter element sealed inside a polypropylene housing. These offer highest packing density of the membrane resulting in a very compact capsule offering long service life.

**Radiation Sterilizable:** AseptiCap KL/KS-

**Autoclavable:** AseptiCap KL/KS




## Specifications


| Construction                                      |   |  |   |                      |
|---|---|--|---|----------------------|
| Final Filter Pore Size                            | 0.1 µm  | 0.2 µm   | 0.45 µm   |                      |
| Pre-filter Pore Size<br>(in case of AseptiCap KS) | 0.45 µm   | 0.8 µm, 0.65 µm, 0.45 µm   |   | 0.8 µm, 0.65 µm      |
| Membrane  | Hydrophilic PES   |  |   |                      |
| Support Layers                                    | Polyester   |  |   |                      |
| Body and Core                                     | Polypropylene   |  |   |                      |
| Integrity Testing/Retention                       |   |  |   |                      |
| Bubble Point                                      | ≥ 31 psi (2.18 Kg/cm <sup>2</sup> )<br>with 50% IPA/Water Solution              | ≥ 50 psi (3.51 Kg/cm <sup>2</sup> )<br>with Water  | ≥ 30 psi (2.11 Kg/cm <sup>2</sup> )<br>with Water                         |                      |
| Microbial Retention                               | LRV >7 for <i>Acholeplasma laidlawii</i><br>(ATCC 23206) per cm <sup>2</sup>    | LRV >7 for <i>Brevundimonas diminuta</i><br>(ATCC 19146) per cm <sup>2</sup>                     | LRV >7 for <i>Serratia marcescens</i><br>(ATCC 14756) per cm <sup>2</sup> |                      |
| Size  |   |  |   |                      |
| Size  | 1"  | 2"   | 5"  | 8"                   |
| Effective Filtration Area (Nominal)               | 250 cm <sup>2</sup>   | 500 cm <sup>2</sup>  | 1000 cm <sup>2</sup>  | 2000 cm <sup>2</sup> |
| Vent and Drain                                    | 1/4" Hose Barb with double platinum cured Silicone 'O' rings for 2", 5" and 8"  |  |   |                      |
| Operational                                       |   |  |   |                      |
| Max. Operating Temperature                        | 0.1 µm  | 80 °C @ ≤ 15 psi (1.0 Kg/cm <sup>2</sup> )   |   |                      |
|   | 0.2 µm  | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )   |   |                      |
| Max. Differential Pressure                        | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C  |  |   |                      |
| Sterilization                                     | By Irradiation  | <b>AseptiCap KL/KS -</b> Gamma Irradiatable up to 50 kGy   |   |                      |
|   | By Gas  | <b>AseptiCap KL/KS:</b> Sterilizable by Ethylene Oxide   |   |                      |
|   |   | <b>AseptiCap KL/KS - :</b> Autoclavable at 121°C for 30 minutes, 1 cycle after gamma irradiation |   |                      |
|   | By Autoclave  | <b>AseptiCap KL/KS:</b> Autoclavable at 121°C for 30 minutes, 25 cycles                          |   |                      |
|   | <b>These cannot be inline steam sterilized</b>                                  |  |   |                      |
| Shelf Life  | 2 years after gamma sterilization<br>3 years after Ethylene Oxide sterilization |  |   |                      |
| pH Compatibility                                  | Compatible with pH range of 1-10  |  |   |                      |

## Water Flow Rates


0.2 µm AseptiCap KS- , 1" Capsule Filters




0.2 µm AseptiCap KS- , 2" Capsule Filters



0.2 µm AseptiCap KS- , 5" Capsule Filters



0.2 µm AseptiCap KS- , 8" Capsule Filters



### End Connection Type:

**A:** ¼" Stepped Hose Barb    **D:** ½" Hose Barb    **E:** 1½" Sanitary Flange    **J:** Quick Connector    **S:** ¾" Sanitary Flange

### Ordering Information

AseptiCap KL/KS and AseptiCap KL/KS -

| Type                            |             | Size |           | Pore Size |           | Inlet/Outlet             |           | Radiation Sterilizable |          | Bell    |          | Sterility     |          | Pack Size |           |
|---------------------------------|-------------|------|-----------|-----------|-----------|--------------------------|-----------|------------------------|----------|---------|----------|---------------|----------|-----------|-----------|
|                                 | Code        |      | Code      |           | Code      |                          | Code      |                        | Code     |         | Code     |               | Code     |           | Code      |
| AseptiCap KL (Single Layer)     | DKLX        | 1"   | 51        | 0.1 µm*   | 36        | ¼" SHB                   | A         | Yes                    | R        | Yes     | B        | Non Sterile   | 1        | 1         | 01        |
| AseptiCap KS (0.45 µm Upstream) | DKSX        | 2"   | 52        | 0.2 µm    | 01        | ¼" MNPT                  | B         | No                     | X        | No Bell | X        | EO Sterile    | 2        |           |           |
| AseptiCap KS (0.65 µm Upstream) | DKS3        | 5"   | 53        | 0.45 µm   | 02        | ½" MNPT                  | C         |                        |          |         |          | Gamma Sterile | 3        |           |           |
| AseptiCap KS (0.8 µm Upstream)  | DKS5        | 8"   | 57        |           |           | ½" Hose Barb             | D         |                        |          |         |          |               |          |           |           |
|                                 |             |      |           |           |           | 1½" Sanitary Flange      | E         |                        |          |         |          |               |          |           |           |
|                                 |             |      |           |           |           | ¾" Sanitary Flange       | S         |                        |          |         |          |               |          |           |           |
|                                 |             |      |           |           |           | Quick Connector          | J         |                        |          |         |          |               |          |           |           |
|                                 |             |      |           |           |           | Single Step ½" Hose Barb | Q         |                        |          |         |          |               |          |           |           |
| <b>Example:</b>                 |             |      |           |           |           |                          |           |                        |          |         |          |               |          |           |           |
|                                 | <b>DKSX</b> |      | <b>57</b> |           | <b>36</b> |                          | <b>DD</b> |                        | <b>R</b> |         | <b>X</b> |               | <b>1</b> |           | <b>01</b> |

\* Note: 0.1 µm is available in AseptiCap KS and AseptiCap KS- only with 0.45 µm upstream

For End Connection, bell availability and dimensions with different sizes refer Page no. 59.

Note: Gamma Sterile filters can not be gamma irradiated again

Example for Non Sterile: DKSX5136DDRDX104

Example for gamma Sterile: DKSX5136DDXX304



## AseptiCap KS- 10", 20", 30"

These are large disposable Polyethersulfone membrane capsule filters for high value biopharma manufacturing processes, providing a unique combination of high throughputs and low hold up volumes. These capsule filters offer serial filtration incorporating a large pore size upstream membrane to protect the downstream membrane and do away with the time and expense associated with assembling, cleaning and validating stainless steel housings.

**Radiation Sterilizable:** AseptiCap KS-

**Autoclavable:** AseptiCap KS




|   |
|---|
| <b>Microbiologically Validated<br/>as per ASTM F 838-05</b> |
| <b>Complies with<br/>USFDA 21 CFR 211.72</b>                |
| <b>Meets and Exceeds<br/>USFDA 21 CFR 177.1520</b>          |

### Specifications


| Construction  |   |  |   |
|---|---|--|---|
| Final Filter Pore Size                              | 0.1 µm  | 0.2 µm   | 0.45 µm   |
| Prefilter Pore Size                                 | 0.45 µm   | 0.8 µm, 0.65 µm, 0.45 µm   | 0.8 µm, 0.65 µm   |
| Membrane  | Hydrophilic PES   |  |   |
| Support Layers                                      | Polyester   |  |   |
| Body and Core                                       | Polypropylene   |  |   |
| Integrity Testing/Retention                         |   |  |   |
| Bubble Point  | ≥ 31 psi (2.18 Kg/cm <sup>2</sup> )<br>with 50% IPA/Water Solution              | ≥ 50 psi (3.51 Kg/cm <sup>2</sup> )<br>with Water  | ≥ 30 psi (2.11 Kg/cm <sup>2</sup> )<br>with Water                         |
| Max. Air Diffusion Flows for<br>10" Capsule Filters | ≤ 29 ml/min @ 50 psi (3.51 Kg/cm <sup>2</sup> )<br>with Water                   | ≤ 30 ml/min @ 37 psi (2.6 Kg/cm <sup>2</sup> )<br>with Water   | ≤ 35 ml/min @ 22 psi (1.54 Kg/cm <sup>2</sup> )<br>with Water             |
| Microbial Retention                                 | LRV >7 for <i>Acholeplasma laidlawii</i><br>(ATCC 23206) per cm <sup>2</sup>    | LRV >7 for <i>Brevundimonas diminuta</i><br>(ATCC 19146) per cm <sup>2</sup>   | LRV >7 for <i>Serratia marcescens</i><br>(ATCC 14756) per cm <sup>2</sup> |
| Size  |   |  |   |
| Size  | 10"   | 20"  | 30"   |
| Effective Filtration Area (Nominal)                 | 6000 cm <sup>2</sup>  | 12000 cm <sup>2</sup>  | 18000 cm <sup>2</sup>   |
| Vent and Drain                                      | 1/4" Hose Barb with double platinum cured Silicone 'O' rings                    |  |   |
| Operational   |   |  |   |
| Max. Operating<br>Temperature                       | 0.1 µm  | 80 °C @ ≤ 15 psi (1.0 Kg/cm <sup>2</sup> )   |   |
|   | 0.2 µm  | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )   |   |
| Max. Differential Pressure                          | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C  |  |   |
| Sterilization                                       | By Irradiation  | <b>AseptiCap KL/KS</b> - Gamma Irradiatable up to 50 kGy   |   |
|   | By Gas  | <b>AseptiCap KL/KS</b> : Sterilizable by Ethylene Oxide  |   |
|   |   | <b>AseptiCap KL/KS</b> - : Autoclavable at 121°C for 30 minutes, 1 cycle after gamma irradiation                           |   |
|   | By Autoclave  | <b>AseptiCap KL/KS</b> : Autoclavable at 121°C for 30 minutes, 25 cycles<br><b>These cannot be inline steam sterilized</b> |   |
| Shelf Life  | 2 years after gamma sterilization<br>3 years after Ethylene Oxide sterilization |  |   |
| pH Compatibility                                    | Compatible with pH range of 1-10  |  |   |

# Water Flow Rates


0.1 µm AseptiCap KS- , 10" Large Capsule Filters




0.2 µm AseptiCap KS- , 10" Large Capsule Filters



0.2 µm AseptiCap KS- , 20" Large Capsule Filters



0.2 µm AseptiCap KS- , 30" Large Capsule Filters



**End Connection Type:** Q: Single Step 1/2" Hose Barb E: 1/2" Sanitary Flange EQ: 1/2" Sanitary Flange Inlet Single Step 1/2" Hose Barb Outlet

## Ordering Information

AseptiCap KS and AseptiCap KS-

| Type                            |      | Size |      | Pore Size |      | Inlet/Outlet               |      | Radiation Sterilizable |      | Inline / T-line |      | Sterility   |      | Pack Size |      |
|---------------------------------|------|------|------|-----------|------|----------------------------|------|------------------------|------|-----------------|------|-------------|------|-----------|------|
|                                 | Code |      | Code |           | Code |                            | Code |                        | Code |                 | Code |             | Code |           | Code |
| AseptiCap KS (0.45 µm Upstream) | LKSX | 10"  | 54   | 0.1 µm*   | 36   | 1/2" Sanitary Flange       | E    | Yes                    | R    | Inline          | X    | Non Sterile | 1    | 1         | 01   |
| AseptiCap KS (0.65 µm Upstream) | LKS3 | 20"  | 55   | 0.2 µm    | 01   | Single Step 1/2" Hose Barb | Q    | No                     | X    | T-line          | T    | EO Sterile  | 2    |           |      |
| AseptiCap KS (0.8 µm Upstream)  | LKS5 | 30"  | 56   | 0.45 µm   | 02   |                            |      | Gamma Sterile          | 3    |                 |      |             |      |           |      |

**Example:**

|      |    |    |    |   |   |   |    |
|------|----|----|----|---|---|---|----|
| LKSX | 54 | 02 | EE | R | T | 1 | 01 |
|------|----|----|----|---|---|---|----|

\* Note: 0.1 µm is available in AseptiCap KS and AseptiCap KS- only with 0.45 µm upstream

For End Connection availability and dimensions with different sizes refer Page no. 59.

Note: Gamma Sterile filters can not be gamma irradiated again

Example for Non Sterile: LKSX5402EERX104

Example for gamma Sterile: LKSX5402EEXX304

# AseptiCap KO- 1", 2", 5", 8"

**AseptiCap KO** capsule filters incorporate a low protein binding PES membrane with polypropylene drainage layers to ensure pH compatibility from 1-14 making these ideal for alkaline fluid streams.



**Microbially Validated  
as per ASTM F 838-05**

**Complies with  
USFDA 21 CFR 211.72**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

## Specifications

| Construction                          |  |   |   |                      |
|---------------------------------------|--|---|---|----------------------|
| Final Filter Pore Size                | 0.2 µm   |   | 0.45 µm   |                      |
| Membrane                              | Hydrophilic PES  |   |   |                      |
| Support Layers                        | Polypropylene  |   |   |                      |
| Body and Core                         | Polypropylene  |   |   |                      |
| Integrity Testing/Retention           |  |   |   |                      |
| Bubble Point                          | ≥ 50 psi (3.51 Kg/cm <sup>2</sup> ) with Water                     |   | ≥ 30 psi (2.11 Kg/cm <sup>2</sup> ) with Water                  |                      |
| Microbial Retention                   | LRV >7 for Brevundimonas diminuta (ATCC 19146) per cm <sup>2</sup> |   | LRV >7 for Serratia marcescens (ATCC 14756) per cm <sup>2</sup> |                      |
| Size                                  |  |   |   |                      |
| Size                                  | 1"   | 2"  | 5"  | 8"                   |
| Effective Filtration Area (Nominal)   | 250 cm <sup>2</sup>  | 500 cm <sup>2</sup>   | 1000 cm <sup>2</sup>  | 2000 cm <sup>2</sup> |
| Vent and Drain                        | 1/4" Hose Barb with double Silicone 'O' rings                      |   |   |                      |
| Operational                           |  |   |   |                      |
| Max. Operating Temperature            | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )                           |   |   |                      |
| Max. Differential Pressure            | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C                             |   |   |                      |
| Sterilization                         | By Gas   | Sterilizable by Ethylene Oxide  |   |                      |
|                                       | By Autoclave   | Autoclavable at 121 °C for 30 minutes, 25 cycles. Can not be in-line steam sterilized |   |                      |
| Typical Water Flow Rates (0.2 µm, 8") | 7.5 lpm @ 0.70 Kg/cm <sup>2</sup> @ 27 °C                          |   |   |                      |
| Shelf Life                            | 3 years after Ethylene Oxide sterilization                         |   |   |                      |
| pH Compatibility                      | Compatible with pH range of 1-14                                   |   |   |                      |

## Ordering Information

| Type            |      | Size |      | Pore Size |      | Inlet/Outlet               |      | X | X | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|----------------------------|------|---|---|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |                            | Code |   |   |             | Code |           | Code |
| AseptiCap KO    | DKLO | 1"   | 51   | 0.2 µm    | 01   | 1/4" SHB                   | A    |   |   | Non Sterile | 1    | 1         | 01   |
|                 |      | 2"   | 52   | 0.45 µm   | 02   | 1/2" Hose Barb             | D    |   |   | EO Sterile  | 2    |           |      |
|                 |      | 5"   | 53   |           |      | 1 1/2" Sanitary Flange     | E    |   |   |             |      |           |      |
|                 |      | 8"   | 57   |           |      | 3/4" Sanitary Flange       | S    |   |   |             |      |           |      |
|                 |      |      |      |           |      | Quick Connector            | J    |   |   |             |      |           |      |
|                 |      |      |      |           |      | Single Step 1/2" Hose Barb | Q    |   |   |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |                            |      |   |   |             |      |           |      |
|                 | DKLO |      | 57   |           | 01   |                            | DD   |   | X |             | X    |           | 1    |
|                 |      |      |      |           |      |                            |      |   |   |             |      |           | 01   |

For End Connection availability and dimensions with different sizes refer Page no. 59.

# AseptiSure HS

**AseptiSure HS** high temperature resistant, serial filtration Polyethersulfone cartridge filters designed to withstand high pressure differential up to 0.3 Kg/cm<sup>2</sup> (5 psi) at high steam sterilization temperatures of up to 135 °C.



**Microbiologically Validated  
as per ASTM F 838-05**

**Complies with  
USFDA 21 CFR 211.72**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

## Specifications

| Construction  |   |  |   |                       |
|---|---|--|---|-----------------------|
| Final Filter Pore Size                              | 0.1 µm  | 0.2 µm   | 0.45 µm   |                       |
| Prefilter Pore Size                                 | 0.45 µm   | 0.8 µm, 0.65 µm, 0.45 µm   | 0.8 µm, 0.65 µm   |                       |
| Membrane  | Hydrophilic PES   |  |   |                       |
| Support Layers                                      | Polyester   |  |   |                       |
| Body and Core                                       | Polypropylene   |  |   |                       |
| Integrity Testing/Retention                         |   |  |   |                       |
| Bubble Point  | ≥ 31 psi (2.18 Kg/cm <sup>2</sup> )<br>with 50% IPA/Water Solution  | ≥ 50 psi (3.51 Kg/cm <sup>2</sup> )<br>with Water                            | ≥ 30 psi (2.11 Kg/cm <sup>2</sup> )<br>with Water                         |                       |
| Air Diffusion Flow<br>(10" Cartridge Filter)        | ≤ 29 ml/min @ 50 psi (3.51 Kg/cm <sup>2</sup> )<br>with Water   | ≤ 30 ml/min @ 37 psi (2.6 Kg/cm <sup>2</sup> )<br>with Water                 | ≤ 35 ml/min @ 22 psi (1.54 Kg/cm <sup>2</sup> )<br>with Water             |                       |
| Microbial Retention                                 | LRV >7 for <i>Acholeplasma laidlawii</i><br>(ATCC 23206) per cm <sup>2</sup>  | LRV >7 for <i>Brevundimonas diminuta</i><br>(ATCC 19146) per cm <sup>2</sup> | LRV >7 for <i>Serratia marcescens</i><br>(ATCC 14756) per cm <sup>2</sup> |                       |
| Size  |   |  |   |                       |
| Size  | 5"  | 10"  | 20"   | 30"                   |
| Effective Filtration Area (Nominal)                 | 3000 cm <sup>2</sup>  | 6000 cm <sup>2</sup>   | 12000 cm <sup>2</sup>   | 18000 cm <sup>2</sup> |
| Operational   |   |  |   |                       |
| Max. Operating Temperature                          | 80 °C @ ≤ 30psi (2 Kg/cm <sup>2</sup> )   |  |   |                       |
| Max. Differential Pressure                          | 50psi (3.5 Kg/cm <sup>2</sup> ) @ 25 °C   |  |   |                       |
| Reverse Pressure                                    | ≤ 10 psi (0.7 Kg/cm <sup>2</sup> ) @ 25 °C  |  |   |                       |
| Sterilization                                       | Inline steam sterilizable upto 135 °C for 30 minutes at a maximum differential pressure of 5 psi (0.3 Kg/cm <sup>2</sup> ), 25 cycles |  |   |                       |
| Typical Water Flow Rates<br>(10" Cartridge filters) | 15 lpm @ 0.70 Kg/cm <sup>2</sup> @ 27 °C  | 24 lpm @ 0.70 Kg/cm <sup>2</sup> @ 27 °C                                     | 32 lpm @ 0.70 Kg/cm <sup>2</sup> @ 27 °C                                  |                       |

## Ordering Information

| Type                                |             | Size |           | Pore Size |           | Adapter |           | Elastomer                    |           | Sterility   |          | Pack Size |           |
|-------------------------------------|-------------|------|-----------|-----------|-----------|---------|-----------|------------------------------|-----------|-------------|----------|-----------|-----------|
|                                     | Code        |      | Code      |           | Code      |         | Code      |                              | Code      |             | Code     |           | Code      |
| AseptiSure HS<br>(0.45 µm Upstream) | CPHX        | 5"   | 53        | 0.1 µm    | 36        | 7P      | A0        | Silicone                     | SS        | Non Sterile | 1        | 1         | 01        |
| AseptiSure HS<br>(0.65 µm Upstream) | CPH3        | 10"  | 54        | 0.2 µm    | 01        | 28      | C0        | Viton                        | SV        |             |          |           |           |
| AseptiSure HS<br>(0.8 µm Upstream)  | CPH5        | 20"  | 55        | 0.45 µm   | 02        | 'O'     | D0        | EPDM                         | SE        |             |          |           |           |
|                                     |             | 30"  | 56        |           |           | 4463    | E0        | FEP<br>Encapsulated<br>Viton | FV*       |             |          |           |           |
|                                     |             |      |           |           |           | 4463B   | H0        |                              |           |             |          |           |           |
| <b>Example:</b>                     |             |      |           |           |           |         |           |                              |           |             |          |           |           |
|                                     | <b>CPHX</b> |      | <b>53</b> |           | <b>01</b> |         | <b>E0</b> |                              | <b>SS</b> |             | <b>1</b> |           | <b>01</b> |

\*FV is available in Adapter Code A0 (7P) only

# AseptiSure KS

**AseptiSure KS** are serial filtration polyethersulfone cartridge filters incorporate a large pore size upstream membrane layer to protect the downstream terminal filtration membrane layer.



**Microbiologically Validated  
as per ASTM F 838-05**

**Complies with  
USFDA 21 CFR 211.72**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

## Specifications

| Construction  |  |  |   |                       |
|---|--|--|---|-----------------------|
| Final Filter Pore Size                              | 0.1 µm   | 0.2 µm   | 0.45 µm   |                       |
| Prefilter Pore Size                                 | 0.45 µm  | 0.8 µm, 0.65 µm, 0.45 µm   | 0.8 µm, 0.65 µm   |                       |
| Membrane  | Hydrophilic PES  |  |   |                       |
| Support Layers                                      | Polyester  |  |   |                       |
| Body and Core                                       | Polypropylene  |  |   |                       |
| Integrity Testing/Retention                         |  |  |   |                       |
| Bubble Point  | ≥ 31 psi (2.18 Kg/cm <sup>2</sup> )<br>with 50% IPA/Water Solution           | ≥ 50 psi (3.51 Kg/cm <sup>2</sup> )<br>with Water                            | ≥ 30 psi (2.11 Kg/cm <sup>2</sup> )<br>with Water                         |                       |
| Air Diffusion Flow<br>(10" Cartridge Filter)        | ≤ 29 ml/min @ 50 psi (3.51 Kg/cm <sup>2</sup> )<br>with Water                | ≤ 30 ml/min @ 37 psi (2.6 Kg/cm <sup>2</sup> )<br>with Water                 | ≤ 35 ml/min @ 22 psi (1.54 Kg/cm <sup>2</sup> )<br>with Water             |                       |
| Microbial Retention                                 | LRV >7 for <i>Acholeplasma laidlawii</i><br>(ATCC 23206) per cm <sup>2</sup> | LRV >7 for <i>Brevundimonas diminuta</i><br>(ATCC 19146) per cm <sup>2</sup> | LRV >7 for <i>Serratia marcescens</i><br>(ATCC 14756) per cm <sup>2</sup> |                       |
| Size  |  |  |   |                       |
| Size  | 5"   | 10"  | 20"   | 30"                   |
| Effective Filtration Area (Nominal)                 | 3000 cm <sup>2</sup>   | 6000 cm <sup>2</sup>   | 12000 cm <sup>2</sup>   | 18000 cm <sup>2</sup> |
| Operational   |  |  |   |                       |
| Max. Operating Temperature                          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )                                     |  |   |                       |
| Max. Differential Pressure                          | 50 psi (3.5 Kg/cm <sup>2</sup> ) @ 25 °C                                     |  |   |                       |
| Reverse Pressure                                    | ≤ 10 psi (0.7 Kg/cm <sup>2</sup> ) @ 25 °C                                   |  |   |                       |
| Sterilization                                       | Autoclavable/In-line steam sterilizable at 121 °C for 30 minutes, 25 cycles  |  |   |                       |
| Typical Water Flow Rates<br>(10" Cartridge filters) | 15 lpm @ 0.70 Kg/cm <sup>2</sup> @ 27 °C                                     | 24 lpm @ 0.70 Kg/cm <sup>2</sup> @ 27 °C                                     | 32 lpm @ 0.70 Kg/cm <sup>2</sup> @ 27 °C                                  |                       |

## Ordering Information

| Type                                |             | Size |           | Pore Size |           | Adapter |           | Elastomer                    |           | Sterility   |          | Pack Size |           |
|-------------------------------------|-------------|------|-----------|-----------|-----------|---------|-----------|------------------------------|-----------|-------------|----------|-----------|-----------|
|                                     | Code        |      | Code      |           | Code      |         | Code      |                              | Code      |             | Code     |           | Code      |
| AseptiSure KS<br>(0.45 µm Upstream) | CPKX        | 5"   | 53        | 0.1 µm    | 36        | 7P      | A0        | Silicone                     | SS        | Non Sterile | 1        | 1         | 01        |
| AseptiSure KS<br>(0.65 µm Upstream) | CPK3        | 10"  | 54        | 0.2 µm    | 01        | 28      | C0        | Viton                        | SV        |             |          |           |           |
| AseptiSure KS<br>(0.8 µm Upstream)  | CPK5        | 20"  | 55        | 0.45 µm   | 02        | 'O'     | D0        | EPDM                         | SE        |             |          |           |           |
|                                     |             | 30"  | 56        |           |           | 4463    | E0        | FEP<br>Encapsulated<br>Viton | FV*       |             |          |           |           |
|                                     |             |      |           |           |           | 4463B   | H0        |                              |           |             |          |           |           |
| <b>Example:</b>                     |             |      |           |           |           |         |           |                              |           |             |          |           |           |
|                                     | <b>CPKX</b> |      | <b>53</b> |           | <b>01</b> |         | <b>E0</b> |                              | <b>SS</b> |             | <b>1</b> |           | <b>01</b> |

\*FV is available in Adapter Code A0 (7P) only

# AseptiSure KR

**AseptiSure KR** cartridge filters incorporate a low protein binding PES membrane with polypropylene drainage layers to ensure pH compatibility from 1-14 making these ideal for alkaline fluid streams.



|   |
|---|
| <b>Microbiologically Validated<br/>as per ASTM F 838-05</b> |
| <b>Complies with<br/>USFDA 21 CFR 211.72</b>                |
| <b>Meets and Exceeds<br/>USFDA 21 CFR 177.1520</b>          |

## Specifications

| Construction                                     |   |  |   |
|--|---|--|---|
| Pore Size  | 0.2 µm  | 0.45 µm  |   |
| Membrane   | Hydrophilic PES   |  |   |
| Support Layers                                   | Polypropylene   |  |   |
| Body and Core                                    | Polypropylene   |  |   |
| Integrity Testing/Retention                      |   |  |   |
| Bubble Point (with Water)                        | ≥ 50 psi (3.51 Kg/cm <sup>2</sup> )   | ≥ 30 psi (2.11 Kg/cm <sup>2</sup> )                                    |   |
| Air Diffusion Flow (10" Cartridge Filter)        | ≤ 40 ml/min @ 37 psi (2.6 Kg/cm <sup>2</sup> ) with Water                   | ≤ 35 ml/min @ 22 psi (1.54 Kg/cm <sup>2</sup> ) with Water             |   |
| Microbial Retention                              | LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm <sup>2</sup>   | LRV >7 for <i>Serratia marcescens</i> (ATCC 14756) per cm <sup>2</sup> |   |
| Size   |   |  |   |
| Size   | 5"  | 10"  | 20" 30"                                     |
| Effective Filtration Area (Nominal)              | 3000 cm <sup>2</sup>  | 6000 cm <sup>2</sup>   | 12000 cm <sup>2</sup> 18000 cm <sup>2</sup> |
| Operational                                      |   |  |   |
| Max. Operating Temperature                       | 80 °C @ ≤ 2 Kg/cm <sup>2</sup> (30 psi)                                     |  |   |
| Max. Differential Pressure                       | 50 psi (3.5 Kg/cm <sup>2</sup> ) @ 25 °C                                    |  |   |
| Reverse Pressure                                 | ≤ 10 psi (0.7 Kg/cm <sup>2</sup> ) @ 25 °C                                  |  |   |
| Sterilization                                    | Autoclavable/In-line steam sterilizable at 121 °C for 30 minutes, 25 cycles |  |   |
| Typical Water Flow Rates (10" Cartridge filters) | 24 lpm @ 0.70 Kg/cm <sup>2</sup> @ 27 °C                                    |  | 32 lpm @ 0.70 Kg/cm <sup>2</sup> @ 27 °C    |

## Ordering Information

| Type            |      | Size |      | Pore Size |      | Adapter |      | Elastomer              |      | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|---------|------|------------------------|------|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |         | Code |                        | Code |             | Code |           | Code |
| AseptiSure KR   | CPKR | 5"   | 53   | 0.2 µm    | 01   | 7P      | A0   | Silicone               | SS   | Non Sterile | 1    | 1         | 01   |
|                 |      | 10"  | 54   | 0.45 µm   | 02   | 28      | C0   | Viton                  | SV   |             |      |           |      |
|                 |      | 20"  | 55   |           |      | 'O'     | D0   | EPDM                   | SE   |             |      |           |      |
|                 |      | 30"  | 56   |           |      | 4463    | E0   | FEP Encapsulated Viton | FV*  |             |      |           |      |
|                 |      |      |      |           |      | 4463B   | H0   |                        |      |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |         |      |                        |      |             |      |           |      |
|                 | CPKR |      | 53   |           | 01   |         | E0   |                        | SS   |             | 1    |           | 01   |

\*FV is available in Adapter Code A0 (7P) only

# Filters for Sterilization of Liquids: Nylon- 66 Membrane Filters

**mdi** Nylon- 66 membrane filters are sterilizing grade filters offering absolute retention and wide chemical compatibility.

**mdi** Nylon filter devices are available as:

| Filter Type           | Single Layer        | Multiple Layer       |
|-----------------------|---------------------|----------------------|
| Capsule Filters       | <i>AseptiCap NL</i> | <i>AseptiCap NS</i>  |
| Cartridge Filters     | –                   | <i>AseptiSure NS</i> |
| Membrane Disc Filters | –                   | NN                   |

## Applications

- Sterilization of compatible solvents and chemicals
- Sterilization of disinfectants in pharmaceutical process and lab areas
- Filtration of hospital disinfectants
- Filtration of rinse water for endoscopes and other hospital equipment and surfaces
- Sterilizing filtration of pharmaceutical, aqueous and non aqueous solutions

## Quality Assurance

These filter devices are manufactured in class 10,000 clean rooms under ISO 9001 : 2008 certified quality management systems and are validated to meet compendia and regulatory requirements.

| Assurance                          |  |
|------------------------------------|--|
| <b>Toxicity</b>                    | Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics  |
| <b>Bioburden</b>                   | Bioburden level is < 1000 cfu/filter device as per ANSI/AAMI/ISO 11737-1 : 1995  |
| <b>Bacterial Endotoxin</b>         | Aqueous extracts exhibit < 0.25 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test  |
| <b>Non Fiber Releasing</b>         | Passes test as per USP and comply with USFDA 21 CFR Part 211.72 for fiber release  |
| <b>Extractables with WFI</b>       | Passes test as per USP   |
| <b>Oxidizable Substances</b>       | Within limits as specified in USP  |
| <b>Particle Shedding</b>           | Passes USP test for particulates in injectables  |
| <b>TOC/Conductivity at 25 °C</b>   | Meets the WFI requirements of USP <643> for Total Organic Carbon and USP <645> for Water Conductivity after a specified volume of purified water flush |
| <b>Indirect Food Additive</b>      | All Polypropylene components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177.1520   |
| <b>Good Manufacturing Practice</b> | These products are manufactured in a facility which adheres to Good Manufacturing Practices  |

## Filter Selection Chart

| Application Area  | Key Application Requirements  | Capsule Filters  | Cartridge Filters   | Disc Filters  |
|---|---|--|---|---|
| Sterilization of compatible solvents and chemicals                          | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- Wide chemical compatibility</li> </ul> | <b>AseptiCap NL/NS</b><br>0.2 µm Nylon- 66 Membrane Capsule Filter | <b>AseptSure NS</b><br>0.2 µm Nylon- 66 Membrane Cartridge Filter | 0.2 µm Nylon- 66 Membrane Disc Filters<br><b>Type- NN</b> |
| Sterilization of disinfectants in pharmaceutical labs and process areas     | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- Wide chemical compatibility</li> </ul> | <b>AseptiCap NL/NS</b><br>0.2 µm Nylon- 66 Membrane Capsule Filter | <b>AseptSure NS</b><br>0.2 µm Nylon-66 Membrane Cartridge Filter  | 0.2 µm Nylon- 66 Membrane Disc Filters<br><b>Type- NN</b> |
| Filtration of hospital disinfectants  | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- Wide chemical compatibility</li> </ul> | <b>AseptiCap NL/NS</b><br>0.2 µm Nylon- 66 Membrane Capsule Filter | <b>AseptSure NS</b><br>0.2 µm Nylon- 66 Membrane Cartridge Filter | 0.2 µm Nylon- 66 Membrane Disc Filters<br><b>Type- NN</b> |
| Filtration of rinse water for endoscopes and other hospital equipments      | <ul style="list-style-type: none"> <li>- Absolute retention</li> </ul>  | <b>AseptiCap NL/NS</b><br>0.2 µm Nylon- 66 Membrane Capsule Filter | <b>AseptSure NS</b><br>0.2 µm Nylon- 66 Membrane Cartridge Filter | -   |
| Sterilizing filtration of pharmaceutical, aqueous and non aqueous solutions | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- Wide chemical compatibility</li> </ul> | <b>AseptiCap NL/NS</b><br>0.2 µm Nylon- 66 Membrane Capsule Filter | <b>AseptSure NS</b><br>0.2 µm Nylon- 66 Membrane Cartridge Filter | 0.2 µm Nylon- 66 Membrane Disc Filters<br><b>Type- NN</b> |



# AseptiCap NL/NS- 1", 2", 5", 8"

**Microbiologically Validated**  
as per ASTM F 838-05

**Complies with**  
**USFDA 21 CFR 211.72**

**Meets and Exceeds**  
**USFDA 21 CFR 177.1520**

**AseptiCap NL/NS-** Nylon- 66 membranes capsule filters are sterilizing grade filters offering absolute retention, wide chemical compatibility, and very low hold up volumes.

**Single Layered:** AseptiCap NL  
**Multiple Layered:** AseptiCap NS


## Specifications

| Construction                                     |  |   |  |
|--|--|---|--|
| Final Filter Pore Size                           | 0.2 µm   | 0.45 µm   |  |
| Prefilter Pore Size<br>(In case of AseptiCap NS) | 0.8 µm, 0.45 µm  | 0.8 µm  |  |
| Membrane   | Nylon- 66  |   |  |
| Support Layers                                   | Polyester  |   |  |
| Body and Core                                    | Polypropylene  |   |  |
| Integrity Testing/Retention                      |  |   |  |
| Bubble Point<br>(with 50% IPA)                   | ≥ 17 psi (1.19 Kg/cm <sup>2</sup> )  | ≥ 11 psi (0.77 Kg/cm <sup>2</sup> )   |  |
| Microbial Retention<br>(LRV >7 for)              | Brevundimonas diminuta<br>(ATCC 19146) per cm <sup>2</sup>                 | Serratia marcescens<br>(ATCC 14756) per cm <sup>2</sup>                       |  |
| Size   |  |   |  |
| Size   | 1"   | 2"  | 5" 8"  |
| Effective Filtration Area<br>(Nominal)           | AseptiCap NL<br>250 cm <sup>2</sup><br>AseptiCap NS<br>200 cm <sup>2</sup> | 900 cm <sup>2</sup><br>700 cm <sup>2</sup>                                    | 1800 cm <sup>2</sup><br>1400 cm <sup>2</sup><br>2700 cm <sup>2</sup><br>2100 cm <sup>2</sup> |
| Vent and Drain                                   | 1/4" Hose Barb with double Silicone 'O' rings                              |   |  |
| Operational                                      |  |   |  |
| Max. Operating Temperature                       | 80 °C @ ≤ 2 Kg/cm <sup>2</sup> (30 psi)                                    |   |  |
| Max. Differential Pressure                       | 4 Kg/cm <sup>2</sup> (60 psi) @ 30 °C                                      |   |  |
| Sterilization                                    | By Gas   | Sterilizable by Ethylene Oxide  |  |
|  | By Autoclave   | Autoclavable at 121 °C for 30 minutes.<br>Can not be in-line steam sterilized |  |



## Water Flow Rates

0.2 µm AseptiCap NL Capsule Filters



## End Connection Type

D: 1/2" Hose Barb

## Ordering Information

| Type                               |      | Size |      | Pore Size |      | Inlet/Outlet               |      | X | Bell    |   | Sterility   |   | Pack Size |    |
|------------------------------------|------|------|------|-----------|------|----------------------------|------|---|---------|---|-------------|---|-----------|----|
|                                    | Code |      | Code |           | Code |                            | Code |   | Code    |   | Code        |   | Code      |    |
| AseptiCap NL                       | DNLX | 1"   | 51   | 0.2 µm    | 01   | 1/4" SHB                   | A    |   | Yes     | B | Non Sterile | 1 | 1         | 01 |
| AseptiCap NS<br>(0.45 µm Upstream) | DNSX | 2"   | 52   | 0.45 µm   | 02   | 1/4" MNPT                  | B    |   | No Bell | X | EO Sterile  | 2 |           |    |
| AseptiCap NS<br>(0.8 µm Upstream)  | DNS5 | 5"   | 53   |           |      | 1/2" MNPT                  | C    |   |         |   |             |   |           |    |
|                                    |      | 8"   | 57   |           |      | 1/2" Hose Barb             | D    |   |         |   |             |   |           |    |
|                                    |      |      |      |           |      | 1 1/2" Sanitary Flange     | E    |   |         |   |             |   |           |    |
|                                    |      |      |      |           |      | 3/4" Sanitary Flange       | S    |   |         |   |             |   |           |    |
|                                    |      |      |      |           |      | Quick Connector            | J    |   |         |   |             |   |           |    |
|                                    |      |      |      |           |      | Single Step 1/2" Hose Barb | Q    |   |         |   |             |   |           |    |
| <b>Example:</b>                    |      |      |      |           |      |                            |      |   |         |   |             |   |           |    |
|                                    | DNSX |      | 52   |           | 01   |                            | DD   |   | X       |   |             | 1 |           | 01 |

For End Connection, bell availability and dimensions with different sizes refer Page no. 59.

# AseptiCap NS- 10", 20", 30"

**Microbiologically Validated**  
as per ASTM F 838-05

**Complies with**  
**USFDA 21 CFR 211.72**

**Meets and Exceeds**  
**USFDA 21 CFR 177.1520**

**AseptiCap NS-** Nylon- 66 membrane large capsule filters are double layered sterilizing grade filters offering absolute retention, wide chemical compatibility, and serial filtration for enhanced throughputs.


## Specifications

| Construction  |   |  |                       |
|---|---|--|-----------------------|
| Final Filter Pore Size                                    | 0.2 µm  | 0.45 µm  |                       |
| Prefilter Pore Size                                       | 0.8 µm, 0.45 µm   | 0.8 µm   |                       |
| Membrane  | Nylon- 66   |  |                       |
| Support Layers  | Polyester   |  |                       |
| Body and Core   | Polypropylene   |  |                       |
| Integrity Testing/Retention                               |   |  |                       |
| Air Diffusion Flow for 10" Capsule Filters (with 50% IPA) | < 25 ml/min @ 16 psi (1.12 Kg/cm <sup>2</sup> )         | < 25 ml/min @ 10 psi (0.7 Kg/cm <sup>2</sup> )                             |                       |
| Microbial Retention (LRV >7 for)                          | Brevundimonas diminuta (ATCC 19146) per cm <sup>2</sup> | Serratia marcescens (ATCC 14756) per cm <sup>2</sup>                       |                       |
| Size  |   |  |                       |
| Size  | 10"   | 20"  | 30"                   |
| Effective Filtration Area (Nominal)                       | 6000 cm <sup>2</sup>                                    | 12000 cm <sup>2</sup>  | 18000 cm <sup>2</sup> |
| Vent and Drain  | 1/4" Hose Barb with double Silicone 'O' rings           |  |                       |
| Operational   |   |  |                       |
| Max. Operating Temperature                                | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )                |  |                       |
| Max. Differential Pressure                                | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C                  |  |                       |
| Sterilization   | By Gas  | Sterilizable by Ethylene Oxide   |                       |
|   | By Autoclave  | Autoclavable at 121 °C for 30 minutes. Can not be in-line steam sterilized |                       |



## Water Flow Rates

0.2 µm AseptiCap NS, 10" Large Capsule Filters



**End Connection Type:** Q: Single Step 1/2" Hose Barb  
E: 1 1/2" Sanitary Flange

## Ordering Information

| Type                            |      | Size |      | Pore Size |      | Inlet/Outlet               |      | X | Inline / T-line |      | Sterility   |      | Pack Size |      |
|---------------------------------|------|------|------|-----------|------|----------------------------|------|---|-----------------|------|-------------|------|-----------|------|
|                                 | Code |      | Code |           | Code |                            | Code |   |                 | Code |             | Code |           | Code |
| AseptiCap NS (0.45 µm Upstream) | LNSX | 10"  | 54   | 0.2 µm    | 01   | 1 1/2" Sanitary Flange     | E    |   | Inline          | X    | Non Sterile | 1    | 1         | 01   |
| AseptiCap NS (0.8 µm Upstream)  | LNS5 | 20"  | 55   | 0.45 µm   | 02   | Single Step 1/2" Hose Barb | Q    |   | T-line          | T    | EO Sterile  | 2    |           |      |
|                                 |      | 30"  | 56   |           |      |                            |      |   |                 |      |             |      |           |      |
| <b>Example:</b>                 |      |      |      |           |      |                            |      |   |                 |      |             |      |           |      |
|                                 | LNSX |      | 54   |           | 01   |                            | QQ   |   |                 | T    |             | 1    |           | 01   |

For End Connection availability and dimensions with different sizes refer Page no. 59.

# AseptiSure NS

**AseptiSure NS**- Nylon- 66 membrane cartridge filters are sterilizing grade filters offering absolute retention and wide chemical compatibility. These filters offer serial filtration for enhanced throughput. The upstream layer is of larger pore size to protect the downstream final filtration layer.

**Microbiologically Validated**  
as per ASTM F 838-05

**Complies with**  
**USFDA 21 CFR 211.72**

**Meets and Exceeds**  
**USFDA 21 CFR 177.1520**




## Specifications

| Construction  |  |  |                       |                       |
|---|--|--|-----------------------|-----------------------|
| Final Filter Pore Size  | 0.2 µm   | 0.45 µm  |                       |                       |
| Prefilter Pore Size   | 0.8 µm, 0.45 µm  | 0.8 µm   |                       |                       |
| Membrane  | Nylon- 66  |  |                       |                       |
| Support Layers  | Polyester  |  |                       |                       |
| Body and Core   | Polypropylene  |  |                       |                       |
| Integrity Testing/Retention                                     |  |  |                       |                       |
| Air Diffusion Flow (with 50% IPA Wetted) (10" Cartridge Filter) | < 25 ml/min @ 16 psi (1.12 Kg/cm <sup>2</sup> )                  | <25 ml/min @ 10 psi (0.7 Kg/cm <sup>2</sup> )        |                       |                       |
| Microbial Retention (LRV >7 for)                                | Brevundimonas diminuta (ATCC 19146) per cm <sup>2</sup>          | Serratia marcescens (ATCC 14756) per cm <sup>2</sup> |                       |                       |
| Size  |  |  |                       |                       |
| Size  | 5"   | 10"  | 20"                   | 30"                   |
| Effective Filtration Area (Nominal)                             | 3000 cm <sup>2</sup>   | 6000 cm <sup>2</sup>                                 | 12000 cm <sup>2</sup> | 18000 cm <sup>2</sup> |
| Operational   |  |  |                       |                       |
| Max. Operating Temperature                                      | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )                         |  |                       |                       |
| Max. Differential Pressure                                      | 50 psi (3.5Kg/cm <sup>2</sup> ) @ 25°C                           |  |                       |                       |
| Reverse Pressure  | < 10 psi (0.7 Kg/cm <sup>2</sup> ) @ 25 °C                       |  |                       |                       |
| Sterilization   | Autoclavable/In-line Steam Sterilizable at 121 °C for 30 minutes |  |                       |                       |

## Water Flow Rates

AseptiSure NS, 10" Cartridge Filters



## Ordering Information

| Type                             |      | Size |      | Pore Size |      | Adapter |      | Elastomer              |      | Sterility   |      | Pack Size |      |
|----------------------------------|------|------|------|-----------|------|---------|------|------------------------|------|-------------|------|-----------|------|
|                                  | Code |      | Code |           | Code |         | Code |                        | Code |             | Code |           | Code |
| AseptiSure NS (0.45 µm Upstream) | CPNX | 5"   | 53   | 0.2 µm    | 01   | 7P      | A0   | Silicone               | SS   | Non Sterile | 1    | 1         | 01   |
| AseptiSure NS (0.8 µm Upstream)  | CPN5 | 10"  | 54   | 0.45 µm   | 02   | 28      | C0   | Viton                  | SV   |             |      |           |      |
|                                  |      | 20"  | 55   |           |      | 'O'     | D0   | EPDM                   | SE   |             |      |           |      |
|                                  |      | 30"  | 56   |           |      | 4463    | E0   | FEP Encapsulated Viton | FV*  |             |      |           |      |
|                                  |      |      |      |           |      | 4463B   | H0   |                        |      |             |      |           |      |

|                 |      |    |    |    |    |   |    |   |   |   |   |   |   |
|-----------------|------|----|----|----|----|---|----|---|---|---|---|---|---|
| <b>Example:</b> | ↓    | ↓  | ↓  | ↓  | ↓  | ↓ | ↓  | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ |
|                 | CPNX | 54 | 01 | A0 | SS | 1 | 01 |   |   |   |   |   |   |

\*FV is available in Adapter Code A0 (7P) only

# Nylon-66 Membrane Disc Filters - Type NN

Microbiologically Validated  
as per ASTM F 838-05

Complies with  
USFDA 21 CFR 211.72

**NN** membrane disc filters are double layered, hydrophilic, non-media migrating, biologically inert, plain white absolute membrane filters offering wide chemical compatibility.



## Specifications

| Construction   |                                       |  |                            |                            |
|--|---------------------------------------|--|----------------------------|----------------------------|
| Pore Size  | 0.2 µm, 0.45 µm, 0.8 µm, 1.2 µm       |  |                            |                            |
| Membrane   | Nylon- 66 polymer                     |  |                            |                            |
| Size   | 90 mm, 142 mm, 293 mm                 |  |                            |                            |
| Integrity Testing/Retention                            |                                       |  |                            |                            |
| Bubble Point<br>(with Water)                           | 0.2 µm                                | ≥ 50 psi (3.51 Kg/cm <sup>2</sup> ) with water         |                            |                            |
|  | 0.45 µm                               | ≥ 32 psi (2.25 Kg/cm <sup>2</sup> ) with water         |                            |                            |
| Microbial<br>Retention                                 | 0.2 µm                                | LRV> 7 for Brevundimonas diminuta as per ASTM F 838-05 |                            |                            |
|  | 0.45 µm                               | LRV> 7 for Serratia marcescens as per ASTM F 838-05    |                            |                            |
| Operational  |                                       |  |                            |                            |
| Max. Operating Temperature                             | 80 °C continuous                      |  |                            |                            |
| Max. Differential Pressure                             | 5 Kg/cm <sup>2</sup>                  |  |                            |                            |
| Water Flow Rates<br>(27 °C @ 0.70 Kg/cm <sup>2</sup> ) | 0.2 µm                                | 0.45 µm  | 0.8 µm                     | 1.2 µm                     |
|  | 14 ml/min/cm <sup>2</sup>             | 37 ml/min/cm <sup>2</sup>                              | 120 ml/min/cm <sup>2</sup> | 180 ml/min/cm <sup>2</sup> |
| Sterilization  | Autoclavable at 121 °C for 30 minutes |  |                            |                            |

## Ordering Information

| Type            |             | Size   |           | Pore Size |           | XX | XX        | Sterility   |           | Pack Size |          |           |
|-----------------|-------------|--------|-----------|-----------|-----------|----|-----------|-------------|-----------|-----------|----------|-----------|
|                 | Code        |        | Code      |           | Code      |    |           |             | Code      |           | Code     |           |
| NN              | NNXX        | 90 mm  | 14        | 0.2 µm    | 01        |    |           | Non Sterile | 1         | 50        | 03       |           |
|                 |             | 142 mm | 16        | 0.45 µm   | 02        |    |           |             |           |           |          |           |
|                 |             | 293 mm | 19        | 0.8 µm    | 03        |    |           |             |           |           |          |           |
|                 |             |        |           | 1.2 µm    | 10        |    |           |             |           |           |          |           |
| <b>Example:</b> |             |        |           |           |           |    |           |             |           |           |          |           |
|                 | <b>NNXX</b> |        | <b>14</b> |           | <b>01</b> |    | <b>XX</b> |             | <b>XX</b> |           | <b>1</b> | <b>03</b> |

# Filters for Sterilization of Air/Gases

## Single Use PVDF Membrane Filter Devices

**mdi** gamma sterilizable *AseptiVent VF* - are hydrophobic PVDF membrane single use capsule filters with a wide range of end connections and different sizes for linear scalability for use with disposable single use assemblies for biopharmaceutical processes.

These filters are validated for microbial retention with liquid bacterial challenge test to ensure reliable performance under worst case conditions.

## Multiple Use PTFE Membrane Filter Devices

**mdi** produces a wide range of PTFE membrane capsule and cartridge filters to meet filtration requirements of biopharmaceutical and pharmaceutical processing.

These filters are validated for microbial retention with liquid bacterial challenge test as per ASTM F838-05 to provide a high degree of sterility assurance for critical applications such as bioreactor/fermentor venting etc. As they offer wide chemical compatibility with organic solvents, these are ideal for manufacture of sterile API.

## Quality Assurance

These filter devices are manufactured in class 10,000 clean rooms under ISO 9001 : 2008 certified quality management systems and are validated to meet compendia and regulatory requirements.

| Filter Type   | Single Use             | Multiple Use         |
|---|------------------------|----------------------|
| Gamma Sterilizable Capsule Filters                              | <i>AseptiVent VF</i> - | -                    |
| Autoclavable Capsule Filters                                    | -                      | <i>AseptiVent TF</i> |
| Steam Sterilizable Cartridge Filters                            | -                      | <i>AseptiSure TF</i> |
| High Temperature Resistant Steam Sterilizable Cartridge Filters | -                      | <i>AseptiSure TH</i> |

## Applications

- Sterile air sparging in fermentors and bioreactors
- Sterile venting of cell factories, bioreactors and fermentors
- Fermentor exhaust
- Sterilization of environmental air in isolators
- Venting of sterile collection vessels
- Cleaning sterile surfaces
- WFI tank venting
- Nitrogen blanketing
- Sterile filtration of API and solvents
- Dry powder injectable filling
- Sterile air for dryers and micronizers

| Assurance                          |  |
|------------------------------------|--|
| <b>Toxicity</b>                    | Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics                              |
| <b>Bioburden</b>                   | Bioburden level is < 1000 cfu/filter device as per ANSI/AAMI/ISO 11737-1 : 1995                        |
| <b>Bacterial Endotoxin</b>         | Aqueous extracts exhibit < 0.5 EU/ml as established by Limulus Amebocyte Lysate (LAL) Test             |
| <b>Non Fiber Releasing</b>         | Passes test as per USP and comply with USFDA 21 CFR Part 211.72 and 210.3 (b)(6) for fiber release     |
| <b>Oxidizable Substances</b>       | Within limits as specified in USP  |
| <b>Particle Shedding</b>           | Passes USP test for particulates in injectables  |
| <b>Indirect Food Additive</b>      | All Polypropylene components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177.1520 |
| <b>Good Manufacturing Practice</b> | These products are manufactured in a facility which adheres to Good Manufacturing Practices            |

## Filter Selection Chart

| Application Area                                   | Key Application Requirements   | Gamma Sterilizable Capsule Filters                             | Steam Sterilizable  |   |
|--|--|--|---|---|
|  |  |  | Capsule Filters   | Cartridge Filters   |
| Sterile venting for fermentors and bioreactors     | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- High hydrophobicity</li> <li>- High flow rates</li> </ul> | <b>AseptiVent VF-y</b><br>0.2 µm PVDF Membrane Capsule Filters | <b>AseptiVent TF</b><br>0.2 µm PTFE Membrane Capsule Filter | <b>AseptiSure TF</b><br>0.2 µm PTFE Membrane Cartridge Filter |
| Sterile air sparging in fermentors and bioreactors | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- High hydrophobicity</li> <li>- High flow rates</li> </ul> | <b>AseptiVent VF-y</b><br>0.2 µm PVDF Membrane Capsule Filters | <b>AseptiVent TF</b><br>0.2 µm PTFE Membrane Capsule Filter | <b>AseptiSure TF</b><br>0.2 µm PTFE Membrane Cartridge Filter |
| Sterile air for cell factories                     | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- High hydrophobicity</li> </ul>                            | <b>AseptiVent VF-y</b><br>0.2µm PVDF Membrane Capsule Filters  | <b>AseptiVent TF</b><br>0.2 µm PTFE Membrane Capsule Filter | -   |
| Venting of small bioreactors                       | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- High hydrophobicity</li> </ul>                            | <b>AseptiVent VF-y</b><br>0.2 µm PVDF Membrane Capsule Filters | <b>AseptiVent TF</b><br>0.2 µm PTFE Membrane Capsule Filter | -   |
| Fermentor exhaust                                  | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- High hydrophobicity</li> <li>- High flow rates</li> </ul> | -  | <b>AseptiVent TF</b><br>0.2 µm PTFE Membrane Capsule Filter | <b>AseptiSure TF</b><br>0.2 µm PTFE Membrane Cartridge Filter |
| Venting of sterile collection vessels              | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- High hydrophobicity</li> <li>- High flow rates</li> </ul> | <b>AseptiVent VF-y</b><br>0.2 µm PVDF Membrane Capsule Filters | <b>AseptiVent TF</b><br>0.2 µm PTFE Membrane Capsule Filter | <b>AseptiSure TF</b><br>0.2 µm PTFE Membrane Cartridge Filter |
| Nitrogen blanketing in sterile API                 | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- High flow rates</li> </ul>                                | -  | <b>AseptiVent TF</b><br>0.2 µm PTFE Membrane Capsule Filter | <b>AseptiSure TF</b><br>0.2 µm PTFE Membrane Cartridge Filter |
| Cleaning sterile surfaces                          | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- High flow rates</li> </ul>                                | -  | <b>AseptiVent TF</b><br>0.2 µm PTFE Membrane Capsule Filter | <b>AseptiSure TF</b><br>0.2 µm PTFE Membrane Cartridge Filter |
| Dry powder injectable filling                      | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- High flow rates</li> </ul>                                | -  | <b>AseptiVent TF</b><br>0.2 µm PTFE Membrane Capsule Filter | <b>AseptiSure TF</b><br>0.2 µm PTFE Membrane Cartridge Filter |
| WFI tank venting                                   | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- High hydrophobicity</li> <li>- High flow rates</li> </ul> | -  | -   | <b>AseptiSure TF</b><br>0.2 µm PTFE Membrane Cartridge Filter |
| Sterile filtration of API/Solvents                 | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- High flow rates</li> </ul>                                | -  | -   | <b>AseptiSure TF</b><br>0.2 µm PTFE Membrane Cartridge Filter |
| Sterile air for dryers and micronizers             | <ul style="list-style-type: none"> <li>- Absolute retention</li> <li>- High flow rates</li> </ul>                                | -  | <b>AseptiVent TF</b><br>0.2 µm PTFE Membrane Capsule Filter | <b>AseptiSure TF</b><br>0.2 µm PTFE Membrane Cartridge Filter |

## AseptiVent VF - , 25mm, 37mm, 50mm

0.2 µm **AseptiVent VF-γ**, gamma sterilizable PVDF membrane vent filters are validated for microbial retention with liquid bacterial challenge test as per ASTM F838-05 to provide a high degree of sterility assurance for critical applications such as small bioreactors, sterile tank venting, bottle venting, barrier filter for vacuum pump etc

**Microbially Validated  
as per ASTM F 838-05**

**Complies with  
USFDA 21 CFR 211.72**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**




### Specifications


| Construction                        |   |   |                        |
|-------------------------------------|---|---|------------------------|
| Final Filter Pore Size              | 0.2 µm  |   |                        |
| Membrane                            | Hydrophobic PVDF  |   |                        |
| Support Layers                      | Polyester   |   |                        |
| Body and Core                       | Polypropylene   |   |                        |
| Integrity Testing/Retention         |   |   |                        |
| Bubble Point                        | ≥ 18 psi (1.27 Kg/cm <sup>2</sup> ) with 50% IPA/ Water Solution                |   |                        |
| Microbial Retention                 | LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm <sup>2</sup>       |   |                        |
| Size                                |   |   |                        |
| Size                                | 25 mm   | 37 mm   | 50 mm                  |
| Effective Filtration Area (Nominal) | 5 cm <sup>2</sup>   | 10 cm <sup>2</sup>  | 20 cm <sup>2</sup>     |
| Operational                         |   |   |                        |
| Max. Operating Temperature          | 60 °C   |   |                        |
| Max. Differential Pressure          | 1.5 Kg/cm <sup>2</sup> (21 psi) @ 30° C   |   |                        |
| Burst Pressure                      | > 14 Kg/cm <sup>2</sup>   | > 8 Kg/cm <sup>2</sup>  | > 8 Kg/cm <sup>2</sup> |
| Sterilization                       | By Irradiation  | Gamma Irradiatable up to 50 kGy   |                        |
|                                     | By Gas  | Sterilizable by Ethylene Oxide  |                        |
|                                     | By Autoclave  | Autoclavable at 121 °C for 30 minutes, 1 cycle after gamma irradiation. Can not be in-line steam sterilized |                        |
| Shelf Life                          | 2 years after gamma sterilization<br>3 years after Ethylene Oxide sterilization |   |                        |

## Air Flow Rates

0.2 µm AseptiVent VF- , 37 mm Capsule Filters



0.2 µm AseptiVent VF- , 50 mm Capsule Filters



### End Connection Type:

**B:** ¼" Stepped Hose Barb

**S:** ¾" Sanitary Flange

## Ordering Information

AseptiVent VF - , 25 mm

| Type            |      | Size  |      | Pore Size |      | Inlet/Outlet     |      | Radiation Sterilizable |      | X | Sterility     |      | Pack Size |      |
|-----------------|------|-------|------|-----------|------|------------------|------|------------------------|------|---|---------------|------|-----------|------|
|                 | Code |       | Code |           | Code |                  | Code |                        | Code |   |               | Code |           | Code |
| AseptiVent VF   | IVFX | 25 mm | 06   | 0.2 µm    | 01   | Female Luer Lock | M    | Yes                    | R    |   | Non Sterile   | 1    | 100       | 04   |
|                 |      |       |      |           |      | Male Luer Slip   | N    |                        |      |   | EO Sterile    | 2    |           |      |
|                 |      |       |      |           |      |                  |      |                        |      |   | Gamma Sterile | 3    |           |      |
| <b>Example:</b> |      |       |      |           |      |                  |      |                        |      |   |               |      |           |      |
|                 | IVFX |       | 06   |           | 01   |                  | MN   |                        | R    |   |               | 1    |           | 04   |

AseptiVent VF - , 37 mm and 50 mm

| Type            |      | Size  |      | Pore Size |      | Inlet/Outlet        |      | Radiation Sterilizable |      | X | Sterility     |      | Pack Size |      |
|-----------------|------|-------|------|-----------|------|---------------------|------|------------------------|------|---|---------------|------|-----------|------|
|                 | Code |       | Code |           | Code |                     | Code |                        | Code |   |               | Code |           | Code |
| AseptiVent VF   | IVFX | 37 mm | 08   | 0.2 µm    | 01   | ¼" SHB              | B    | Yes                    | R    |   | Non Sterile   | 1    | 12        | 08   |
|                 |      | 50 mm | 10   |           |      | ¾" Sanitary Flange* | S    |                        |      |   | EO Sterile    | 2    | 20        | 09   |
|                 |      |       |      |           |      |                     |      |                        |      |   | Gamma Sterile | 3    |           |      |
| <b>Example:</b> |      |       |      |           |      |                     |      |                        |      |   |               |      |           |      |
|                 | IVFX |       | 08   |           | 01   |                     | BB   |                        | R    |   |               | 1    |           | 09   |

\* Available only in 50 mm



## AseptiVent VF - , 1", 2", 5", 8"

**Microbially Validated  
as per ASTM F 838-05**

**Complies with  
USFDA 21 CFR 211.72**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

**AseptiVent VF -γ** PVDF membrane vent filters are validated for microbial retention with liquid bacterial challenge test as per ASTM F838-05 to provide a high degree of sterility assurance for critical applications.




### Specifications


| Construction                        |   |  |                      |                      |
|-------------------------------------|---|--|----------------------|----------------------|
| Pore Size                           | 0.2 μm  |  |                      |                      |
| Membrane                            | Hydrophobic PVDF  |  |                      |                      |
| Support Layers                      | Polyester   |  |                      |                      |
| Body and Core                       | Polypropylene   |  |                      |                      |
| Integrity Testing/Retention         |   |  |                      |                      |
| Bubble Point                        | ≥ 18 psi (1.26 Kg/cm <sup>2</sup> ) with 50% IPA/Water Solution                 |  |                      |                      |
| Microbial Retention                 | LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm <sup>2</sup>       |  |                      |                      |
| Size                                |   |  |                      |                      |
| Size                                | 1"  | 2"   | 5"                   | 8"                   |
| Effective Filtration Area (Nominal) | 250 cm <sup>2</sup>   | 500 cm <sup>2</sup>  | 1000 cm <sup>2</sup> | 2000 cm <sup>2</sup> |
| Vent and Drain                      | 1/4" Hose Barb with double Silicone 'O' rings for 2", 5" and 8" capsule filters |  |                      |                      |
| Operational                         |   |  |                      |                      |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )  |  |                      |                      |
| Max. Differential Pressure          | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C  |  |                      |                      |
| Sterilization                       | By Irradiation  | Gamma Irradiatable up to 50 kGy  |                      |                      |
|                                     | By Gas  | Sterilizable by Ethylene Oxide   |                      |                      |
|                                     | By Autoclave  | Autoclavable at 121°C for 30 minutes, 1 cycle after gamma irradiation. Can not be in-line steam sterilized |                      |                      |
| Shelf Life                          | 2 years after gamma sterilization   |  |                      |                      |
|                                     | 3 years after Ethylene Oxide sterilization                                      |  |                      |                      |

## Air Flow Rates


0.2 µm AseptiVent VF- , 1" Capsule Filters, DD



0.2 µm AseptiVent VF- , 2" Capsule Filters, DD



0.2 µm AseptiVent VF- , 5" Capsule Filters, DD



### End Connection Type

D: ½" Hose Barb

## Ordering Information

| Type            |             | Size      |           | Pore Size |      | Inlet/Outlet             |      | Radiation Sterilizable |      | X        | Sterility     |           | Pack Size |      |
|-----------------|-------------|-----------|-----------|-----------|------|--------------------------|------|------------------------|------|----------|---------------|-----------|-----------|------|
|                 | Code        |           | Code      |           | Code |                          | Code |                        | Code |          |               |           | Code      | Code |
| AseptiVent VF   | DVLX        | 1"        | 51        | 0.2 µm    | 01   | ¼" SHB                   | A    | Yes                    | R    |          | Non Sterile   | 1         | 1         | 01   |
|                 |             | 2"        | 52        |           |      | ½" Hose Barb             | D    |                        |      |          | EO Sterile    | 2         |           |      |
|                 |             | 5"        | 53        |           |      | 1½" Sanitary Flange      | E    |                        |      |          | Gamma Sterile | 3         |           |      |
|                 |             | 8"        | 57        |           |      | ¾" Sanitary Flange       | S    |                        |      |          |               |           |           |      |
|                 |             |           |           |           |      | Quick Connector          | J    |                        |      |          |               |           |           |      |
|                 |             |           |           |           |      | Single Step ½" Hose Barb | Q    |                        |      |          |               |           |           |      |
| <b>Example:</b> | <b>DVLX</b> | <b>57</b> | <b>01</b> |           |      | <b>EE</b>                |      | <b>R</b>               |      | <b>X</b> | <b>1</b>      | <b>01</b> |           |      |

For End Connection availability and dimensions with different sizes refer Page no. 59.

## AseptiVent TF- 25 mm, 37 mm, 50 mm

**Microbiologically Validated  
as per ASTM F 838-05**

**Complies with  
USFDA 21 CFR 211.72**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

**AseptiVent TF** Disposable inline PTFE gas filters are convenient pre-fabricated devices used for sterilization of gases and as a bacterial air vent in various pharmaceutical and biopharmaceutical processes.




### Specifications


| Construction                        |   |   |   |
|-------------------------------------|---|---|---|
| Final Filter Pore Size              | 0.2 µm  |   | 0.45 µm   |
| Membrane                            | Hydrophobic PTFE  |   |   |
| Support Layers                      | Polypropylene   |   |   |
| Body and Core                       | Polypropylene   |   |   |
| Integrity Testing/Retention         |   |   |   |
| Bubble Point                        | ≥ 22 psi (1.54 Kg/cm <sup>2</sup> ) with 70% IPA/Water Solution           |   | ≥ 10 psi (0.7 Kg/cm <sup>2</sup> ) with 70% IPA/Water Solution        |
| Microbial Bacterial Retention       | LRV >7 for <i>Brevundimonas diminuta</i> (ATCC 19146) per cm <sup>2</sup> |   | LRV >7 for <i>Serratia marcescens</i> ATCC 14756) per cm <sup>2</sup> |
| Size                                |   |   |   |
| Size                                | 25 mm   | 37 mm   | 50 mm   |
| Effective Filtration Area (Nominal) | 5 cm <sup>2</sup>   | 10 cm <sup>2</sup>  | 20 cm <sup>2</sup>  |
| Operational                         |   |   |   |
| Max. Operating Temperature          | 60 °C   |   |   |
| Max. Differential Pressure          | 42 psi (3 Kg/cm <sup>2</sup> ) @ 30 °C                                    |   |   |
| Burst Pressure                      | > 14 Kg/cm <sup>2</sup>   | > 8 Kg/cm <sup>2</sup>  | > 8 Kg/cm <sup>2</sup>  |
| Sterilization                       | By Gas  | Sterilizable by Ethylene Oxide  |   |
|                                     | By Autoclave  | Autoclavable at 121 °C for 30 minutes, 30 cycles. Can not be in-line steam sterilized |   |
| Shelf Life                          | 3 years after Ethylene Oxide sterilization                                |   |   |

## Air Flow Rates


0.2 µm AseptiVent TF, 25 mm Capsule Filters



0.2 µm AseptiVent TF, 37 mm Capsule Filters



0.2 µm AseptiVent TF, 50 mm Capsule Filters



## End Connection Type:

B: ¼" Stepped Hose Barb C: 1/8" MNPT S: ¾" Sanitary Flange

## Ordering Information

AseptiVent TF- 25 mm

| Type          |      | Size  |      | Pore Size |      | Inlet/Outlet     |      | XX | Sterility   |      | Pack Size |      |
|---------------|------|-------|------|-----------|------|------------------|------|----|-------------|------|-----------|------|
|               | Code |       | Code |           | Code |                  | Code |    |             | Code |           | Code |
| AseptiVent TF | ITFX | 25 mm | 06   | 0.2 µm    | 01   | Female Luer Lock | M    |    | Non Sterile | 1    | 100       | 04   |
|               |      |       |      | 0.45 µm   | 02   | Male Luer Slip   | N    |    | EO Sterile  | 2    |           |      |

**Example:** ITFX 06 01 MN XX 1 04

AseptiVent TF- 37 mm, 50 mm

| Type          |      | Size  |      | Pore Size |      | Inlet/Outlet        |      | XX | Sterility   |      | Pack Size |      |
|---------------|------|-------|------|-----------|------|---------------------|------|----|-------------|------|-----------|------|
|               | Code |       | Code |           | Code |                     | Code |    |             | Code |           | Code |
| AseptiVent TF | ITFX | 37 mm | 08   | 0.2 µm    | 01   | *¼" SHB             | B    |    | Non Sterile | 1    | 12        | 08   |
|               |      | 50 mm | 10   | 0.45 µm   | 02   | 1/8" MNPT           | C    |    | EO Sterile  | 2    | 20        | 09   |
|               |      |       |      |           |      | ¾" Sanitary Flange* | S    |    |             |      |           |      |

**Example:** ITFX 08 01 BB XX 1 09

\* Note: AseptiVent TF- 37 mm is available with BB connection only

# AseptiVent TF- 1", 2", 5", 8"

**Microbiologically Validated**  
as per ASTM F 838-05

**Complies with**  
**USFDA 21 CFR 211.72**

**Meets and Exceeds**  
**USFDA 21 CFR 177.1520**

**AseptiVent TF** capsule filters employ hydrophobic PTFE membrane offering absolute retention and very wide chemical compatibility making these useful for sterile filtration of air/gases as well as aggressive solvents.


## Specifications

| Construction                        |   |   |  |                      |
|-------------------------------------|---|---|--|----------------------|
| Final Filter Pore Size              | 0.2 µm  |   | 0.45 µm  |                      |
| Membrane                            | Hydrophobic PTFE  |   |  |                      |
| Support Layers                      | Polypropylene   |   |  |                      |
| Body and Core                       | Polypropylene   |   |  |                      |
| Integrity Testing/Retention         |   |   |  |                      |
| Bubble Point (with 70% IPA Wetted)  | ≥ 22 psi (1.55 Kg/cm <sup>2</sup> )   |   | ≥ 10 psi (0.7 Kg/cm <sup>2</sup> )                   |                      |
| Microbial Retention (LRV >7 for)    | Brevundimonas diminuta (ATCC 19146) per cm <sup>2</sup>                         |   | Serratia marcescens (ATCC 14756) per cm <sup>2</sup> |                      |
| Size                                |   |   |  |                      |
| Size                                | 1"  | 2"  | 5"   | 8"                   |
| Effective Filtration Area (Nominal) | 250 cm <sup>2</sup>   | 500 cm <sup>2</sup>   | 1000 cm <sup>2</sup>                                 | 2000 cm <sup>2</sup> |
| Vent and Drain                      | 1/4" Hose Barb with double Silicone 'O' rings for 2", 5" and 8" capsule filters |   |  |                      |
| Operational                         |   |   |  |                      |
| Max. Operating Temperature          | 80 °C @ ≤30 psi (2 Kg/cm <sup>2</sup> )   |   |  |                      |
| Max. Differential Pressure          | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C  |   |  |                      |
| Sterilization                       | By Gas  | Sterilizable by Ethylene Oxide  |  |                      |
|                                     | By Autoclave  | Autoclavable at 121 °C for 30 minutes, 50 cycles. Can not be in-line steam sterilized |  |                      |
| Shelf Life                          | 3 years after Ethylene Oxide sterilization                                      |   |  |                      |



## Air Flow Rates

0.2 µm AseptiVent TF Capsule Filters



## End Connection Type

D: ½" Hose Barb

## Ordering Information

| Type            |      | Size |      | Pore Size |      | Inlet/Outlet             |      | X | X | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|--------------------------|------|---|---|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |                          | Code |   |   |             | Code |           | Code |
| AseptiVent TF   | DTLX | 1"   | 51   | 0.2 µm    | 01   | ¼" SHB                   | A    |   |   | Non Sterile | 1    | 1         | 01   |
|                 |      | 2"   | 52   | 0.45 µm   | 02   | ¼" MNPT                  | B    |   |   | EO Sterile  | 2    |           |      |
|                 |      | 5"   | 53   |           |      | ½" MNPT                  | C    |   |   |             |      |           |      |
|                 |      | 8"   | 57   |           |      | ½" Hose Barb             | D    |   |   |             |      |           |      |
|                 |      |      |      |           |      | 1½" Sanitary Flange      | E    |   |   |             |      |           |      |
|                 |      |      |      |           |      | ¾" Sanitary Flange       | S    |   |   |             |      |           |      |
|                 |      |      |      |           |      | Quick Connector          | J    |   |   |             |      |           |      |
|                 |      |      |      |           |      | Single Step ½" Hose Barb | Q    |   |   |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |                          |      |   |   |             |      |           |      |
|                 | DTLX |      | 53   |           | 01   |                          | DD   |   | X |             | X    |           | 1    |
|                 |      |      |      |           |      |                          |      |   |   |             |      |           | 01   |

For End Connection availability and dimensions with different sizes refer Page no. 59.

# AseptiVent TF- 10", 20", 30"

**AseptiVent TF** PTFE large capsule filters offer absolute retention and wide chemical compatibility for sterile filtration of air/gases as well as aggressive solvents in large volume.

**Microbiologically Validated**  
as per ASTM F 838-05

**Complies with**  
**USFDA 21 CFR 211.72**

**Meets and Exceeds**  
**USFDA 21 CFR 177.1520**




## Specifications

| Construction  |   |   |                       |
|---|---|---|-----------------------|
| Pore Size   | 0.2 µm  | 0.45 µm   |                       |
| Membrane  | Hydrophobic PTFE  |   |                       |
| Support Layers  | Polypropylene   |   |                       |
| Body and Core   | Polypropylene   |   |                       |
| Integrity Testing/Retention                               |   |   |                       |
| Air Diffusion Flow (70% IPA Wetted) (10" Capasule Filter) | ≤ 45 ml/min @ 16 psi (1.12 Kg/cm <sup>2</sup> )         | ≤ 45 ml/min @ 8 psi (0.56 Kg/cm <sup>2</sup> )  |                       |
| Microbial Retention (LRV >7 for)                          | Brevundimonas diminuta (ATCC 19146) per cm <sup>2</sup> | Serratia marcescens (ATCC 14756) per cm <sup>2</sup>                                  |                       |
| Size  |   |   |                       |
| Size  | 10"   | 20"   | 30"                   |
| Effective Filtration Area (Nominal)                       | 6000 cm <sup>2</sup>                                    | 12000 cm <sup>2</sup>   | 18000 cm <sup>2</sup> |
| Vent and Drain  | 1/4" Hose Barb with double Silicone 'O' rings           |   |                       |
| Operational   |   |   |                       |
| Max. Operating Temperature                                | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )                |   |                       |
| Max. Differential Pressure                                | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C                  |   |                       |
| Sterilization   | By Gas  | Sterilizable by Ethylene Oxide  |                       |
|   | By Autoclave  | Autoclavable at 121 °C for 30 minutes, 30 cycles. Can not be in-line steam sterilized |                       |
| Shelf Life  | 3 years after Ethylene Oxide sterilization              |   |                       |

## Air Flow Rates

0.2 µm AseptiVent TF, 10" Large Capsule Filters



## End Connection Type

E: 1½" Sanitary Flange

## Ordering Information

| Type            |             | Size      |           | Pore Size |          | Inlet/Outlet             |          | X         | Inline / T-line |      | Sterility   |      | Pack Size |      |
|-----------------|-------------|-----------|-----------|-----------|----------|--------------------------|----------|-----------|-----------------|------|-------------|------|-----------|------|
|                 | Code        |           | Code      |           | Code     |                          | Code     |           |                 | Code |             | Code |           | Code |
| AseptiVent TF   | LTLX        | 10"       | 54        | 0.2 µm    | 01       | 1½" Sanitary Flange      | E        |           | Inline          | X    | Non Sterile | 1    | 1         | 01   |
|                 |             | 20"       | 55        | 0.45 µm   | 02       | Single Step ½" Hose Barb | Q        |           | T-line          | T    | EO Sterile  | 2    |           |      |
|                 |             | 30"       | 56        |           |          |                          |          |           |                 |      |             |      |           |      |
| <b>Example:</b> | <b>LTLX</b> | <b>54</b> | <b>01</b> | <b>EE</b> | <b>X</b> | <b>X</b>                 | <b>1</b> | <b>01</b> |                 |      |             |      |           |      |

For End Connection availability and dimensions with different sizes refer Page no. 59.

# AseptiSure TH

**Microbiologically Validated  
as per ASTM F 838-05**

**Complies with  
USFDA 21 CFR 211.72**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

**AseptiSure TH** cartridge filters are high temperature resistant PTFE filters which are steam sterilizable at upto 135°C. These filters are validated with liquid microbial challenge test as per ASTM F 838-05 to offer absolute retention even under high moisture conditions.


## Specifications

| Construction  |  |  |                      |                       |
|---|--|--|----------------------|-----------------------|
| Final Filter Pore Size  | 0.2 µm   |  | 0.45 µm              |                       |
| Membrane  | Hydrophobic PTFE   |  |                      |                       |
| Support Layers  | Polypropylene  |  |                      |                       |
| Body and Core   | Polypropylene  |  |                      |                       |
| Integrity Testing/Retention                                     |  |  |                      |                       |
| Air Diffusion Flow (with 70% IPA Wetted) (10" Cartridge Filter) | ≤ 45 ml/min @ 16 psi (1.12 Kg/cm <sup>2</sup> )  | ≤ 45 ml/min @ 8 psi (0.56 Kg/cm <sup>2</sup> )       |                      |                       |
| Microbial Retention (LRV >7 for)                                | Brevundimonas diminuta (ATCC 19146) per cm <sup>2</sup>  | Serratia marcescens (ATCC 14756) per cm <sup>2</sup> |                      |                       |
| Size  |  |  |                      |                       |
| Size  | 5"   | 10"  | 20"                  | 30"                   |
| Effective Filtration Area (Nominal)                             | 3000 cm <sup>2</sup>   | 6000 cm <sup>2</sup>                                 | 1200 cm <sup>2</sup> | 18000 cm <sup>2</sup> |
| Operational   |  |  |                      |                       |
| Max. Operating Temperature                                      | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )   |  |                      |                       |
| Max. Differential Pressure                                      | 50 psi (3.5 Kg/cm <sup>2</sup> ) @ 25 °C   |  |                      |                       |
| Reverse Pressure  | ≤ 10 psi (0.7 Kg/cm <sup>2</sup> ) @ 25 °C   |  |                      |                       |
| Sterilization   | Autoclavable/In-line steam sterilizable upto 135 °C, 80 cycles @ maximum differential pressure of 5 psi (0.3 Kg/cm <sup>2</sup> ) for 30minutes. |  |                      |                       |



## Air Flow Rates

0.2 µm AseptiSure TH Cartridge Filters



## Ordering Information

| Type            |      | Size |      | Pore Size |      | Adapter |      | Elastomer              |      | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|---------|------|------------------------|------|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |         | Code |                        | Code |             | Code |           | Code |
| AseptiSure TH   | CPTH | 5"   | 53   | 0.2 µm    | 01   | 7P      | A0   | Silicone               | SS   | Non Sterile | 1    | 1         | 01   |
|                 |      | 10"  | 54   | 0.45 µm   | 02   | 28      | C0   | Viton                  | SV   |             |      |           |      |
|                 |      | 20"  | 55   |           |      | 'O'     | D0   | EPDM                   | SE   |             |      |           |      |
|                 |      | 30"  | 56   |           |      | 4463    | E0   | FEP Encapsulated Viton | FV*  |             |      |           |      |
|                 |      |      |      |           |      | 4463B   | H0   |                        |      |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |         |      |                        |      |             |      |           |      |
|                 | CPTH |      | 56   |           | 01   |         | A0   |                        | SS   |             | 1    |           | 01   |

\*FV is available in Adapter Code A0 (7P) only

# AseptiSure TF

**Microbially Validated  
as per ASTM F 838-05**

**Complies with  
USFDA 21 CFR 211.72**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

**AseptiSure TF** cartridge filters employ hydrophobic PTFE membrane offering absolute retention, wide chemical compatibility, and are validated with liquid bacterial challenge test.




## Specifications

| Construction  |  |  |   |
|---|--|--|---|
| Final Filter Pore Size  | 0.2 µm   | 0.45 µm  |   |
| Membrane  | Hydrophobic PTFE   |  |   |
| Support Layers  | Polypropylene  |  |   |
| Body and Core   | Polypropylene  |  |   |
| Integrity Testing/Retention                                     |  |  |   |
| Air Diffusion Flow (with 70% IPA Wetted) (10" Cartridge Filter) | ≤ 45 ml/min @ 16 psi (1.12 Kg/cm <sup>2</sup> )                              | ≤ 45 ml/min @ 8 psi (0.56 Kg/cm <sup>2</sup> )       |   |
| Microbial Retention (LRV >7 for)                                | Brevundimonas diminuta (ATCC 19146) per cm <sup>2</sup>                      | Serratia marcescens (ATCC 14756) per cm <sup>2</sup> |   |
| Size  |  |  |   |
| Size  | 5"   | 10"  | 20" 30"                                     |
| Effective Filtration Area (Nominal)                             | 3000 cm <sup>2</sup>   | 6000 cm <sup>2</sup>                                 | 12000 cm <sup>2</sup> 18000 cm <sup>2</sup> |
| Operational   |  |  |   |
| Max. Operating Temperature                                      | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )                                     |  |   |
| Max. Differential Pressure                                      | 50 psi (3.5 Kg/cm <sup>2</sup> ) @ 25 °C                                     |  |   |
| Reverse Pressure  | ≤ 10 psi (0.7 Kg/cm <sup>2</sup> ) @ 25 °C                                   |  |   |
| Sterilization   | Autoclavable/In-line steam sterilizable at 121 °C for 30 minutes, 100 cycles |  |   |

## Air Flow Rates

0.2 µm AseptiSure TF Cartridge Filters



## Ordering Information

| Type            |      | Size |      | Pore Size |      | Adapter |      | Elastomer              |      | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|---------|------|------------------------|------|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |         | Code |                        | Code |             | Code |           | Code |
| AseptiSure TF   | CPTF | 5"   | 53   | 0.2 µm    | 01   | 7P      | A0   | Silicone               | SS   | Non Sterile | 1    | 1         | 01   |
|                 |      | 10"  | 54   | 0.45 µm   | 02   | 28      | C0   | Viton                  | SV   |             |      |           |      |
|                 |      | 20"  | 55   |           |      | 'O'     | D0   | EPDM                   | SE   |             |      |           |      |
|                 |      | 30"  | 56   |           |      | 4463    | E0   | FEP Encapsulated Viton | FV*  |             |      |           |      |
|                 |      |      |      |           |      | 4463B   | H0   |                        |      |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |         |      |                        |      |             |      |           |      |
|                 | CPTF |      | 54   |           | 01   |         | A0   |                        | FV   |             | 1    |           | 01   |

\*FV is available in Adapter Code A0 (7P) only



# Filters for Clarification and Pre-filtration

**mdi** offers a range of pre-filters designed to protect the terminal sterilizing grade membrane filters and maximize throughputs.

These are biologically and chemically inert filters offering very high retention efficiency and are available in cartridge filter and capsule filter formats, in different sizes, pore sizes, and end connections to suit different needs.

These filter devices are available as:

| Filter Type   | Product Name  |
|---|---|
| Polyethersulfone Membrane Capsule Filters with Microglassfiber Upstream   | <i>ClariPro GK- 1", 2", 5", 8"</i><br><i>ClariPro GK- 10", 20", 30"</i> |
| Polyethersulfone Membrane Cartridge Filters with Microglassfiber Upstream | <i>ClariSure GK</i>   |
| Microglassfiber Capsule Filters   | <i>ClariCap GS- 1", 2", 5", 8"</i><br><i>ClariCap GS- 10", 20", 30"</i> |
| Microglassfiber Cartridge Filters   | <i>ClariSure GS</i><br><i>ClariSure GP</i>                              |
| Polypropylene Capsule Filters   | <i>ClariCap PP-1", 2", 5", 8"</i><br><i>ClariCap PP-10", 20", 30"</i>   |
| Polypropylene Cartridge Filters   | <i>ClariSure PA</i>   |

## Applications

- Precipitate removal post viral inactivation
- Pre-filtration of cell culture media
- Pre-filtration of serum and other viscous biologicals
- Pre-filtration of serum solutions
- Clarification of cell harvest supernatant
- Pre-filtration of protein solutions
- Pre-filtration of high value difficult to filter drug solutions
- Pre-filtration of large volume parenterals
- Pre-filtration of difficult to filter SVP
- Polishing of turbid solutions
- Pre-filtration of fermentor air

## Quality Assurance

These filter devices are manufactured in class 10,000 clean rooms under ISO 9001 : 2008 certified quality management systems and are deeply validated to meet compendia and regulatory requirements.

| Assurance                          |  |
|------------------------------------|--|
| <b>Toxicity</b>                    | Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics                              |
| <b>Non Fiber Releasing</b>         | Passes test as per USP and comply with USFDA 21 CFR Part 210.3 (b)(6) for fiber release                |
| <b>Extractables with WFI</b>       | Passes test as per USP   |
| <b>Oxidizable Substances</b>       | Within limits as specified in USP  |
| <b>Particle Shedding</b>           | Passes USP test for particulates in injectables  |
| <b>Indirect Food Additive</b>      | All Polypropylene components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177.1520 |
| <b>Good Manufacturing Practice</b> | These products are manufactured in a facility which adheres to Good Manufacturing Practices            |

# Filter Selection Chart

| Application Area  | Key Application Requirements  | Capsule Filters  |   | Cartridge Filters   |  |
|---|---|--|---|---|--|
| <b>Biopharmaceuticals</b>                                       |   |  |   |   |  |
| Precipitate removal post viral inactivation                     | <ul style="list-style-type: none"> <li>- High retention efficiency</li> <li>- High throughput</li> </ul>                              | <b>ClariPro GK</b><br>PES Membrane Capsule Filters with Microglassfiber upstream | <b>ClariCap GS</b><br>Microglassfiber Capsule Filters | <b>ClariSure GK</b><br>PES Membrane Cartridge Filters with Microglassfiber upstream | <b>ClariSure GS</b><br>Microglassfiber Cartridge Filters |
| Pre-filtration of cell culture media                            | <ul style="list-style-type: none"> <li>- High retention efficiency</li> <li>- High throughput</li> </ul>                              | <b>ClariPro GK</b><br>PES Membrane Capsule Filters with Microglassfiber upstream | <b>ClariCap GS</b><br>Microglassfiber Capsule Filters | <b>ClariSure GK</b><br>PES Membrane Cartridge Filters with Microglassfiber upstream | <b>ClariSure GS</b><br>Microglassfiber Cartridge Filters |
| Prefiltration of serum and other viscous biologicals            | <ul style="list-style-type: none"> <li>- High retention efficiency</li> <li>- High throughput</li> </ul>                              | <b>ClariPro GK</b><br>PES Membrane Capsule Filters with Microglassfiber upstream | <b>ClariCap GS</b><br>Microglassfiber Capsule Filters | <b>ClariSure GK</b><br>PES Membrane Cartridge Filters with Microglassfiber upstream | <b>ClariSure GS</b><br>Microglassfiber Cartridge Filters |
| Pre-filtration of serum solutions                               | <ul style="list-style-type: none"> <li>- High retention efficiency</li> <li>- High throughput</li> </ul>                              | <b>ClariPro GK</b><br>PES Membrane Capsule Filters with Microglassfiber upstream | <b>ClariCap GS</b><br>Microglassfiber Capsule Filters | <b>ClariSure GK</b><br>PES Membrane Cartridge Filters with Microglassfiber upstream | <b>ClariSure GS</b><br>Microglassfiber Cartridge Filters |
| Clarification of cell harvest supernatants                      | <ul style="list-style-type: none"> <li>- High retention efficiency</li> <li>- High throughput</li> </ul>                              | <b>ClariPro GK</b><br>PES Membrane Capsule Filters with Microglassfiber upstream | <b>ClariCap GS</b><br>Microglassfiber Capsule Filters | <b>ClariSure GK</b><br>PES Membrane Cartridge Filters with Microglassfiber upstream | <b>ClariSure GS</b><br>Microglassfiber Cartridge Filters |
| Pre-filtration of fermentor air                                 | <ul style="list-style-type: none"> <li>- High retention efficiency</li> </ul>   | -  | <b>ClariCap PP</b><br>Polypropylene Capaule Filters   | <b>ClariSure PA</b><br>Pleated Polypropylene Cartridge Filters                      | -  |
| Pre-filtration of proteinaceous liquids                         | <ul style="list-style-type: none"> <li>- Low hold up volume</li> <li>- High throughput</li> </ul>                                     | <b>ClariPro GK</b><br>PES Membrane Capsule Filters with Microglassfiber upstream | <b>ClariCap PP</b><br>Polypropylene Capaule Filters   | -   | -  |
| <b>Pharmaceuticals</b>  |   |  |   |   |  |
| Pre-filtration of high value difficult to filter drug solutions | <ul style="list-style-type: none"> <li>- High retention efficiency</li> <li>- High throughput</li> </ul>                              | <b>ClariPro GK</b><br>PES Membrane Capsule Filters with Microglassfiber upstream | <b>ClariCap GS</b><br>Microglassfiber Capsule Filters | <b>ClariSure GK</b><br>PES Membrane Cartridge Filters with Microglassfiber upstream | <b>ClariSure GS</b><br>Microglassfiber Cartridge Filters |
| Pre-filtration of large volume parenterals                      | <ul style="list-style-type: none"> <li>- High retention efficiency</li> <li>- High throughput</li> </ul>                              | -  | -   | <b>ClariSure PA</b><br>Pleated Polypropylene Cartridge Filters                      | <b>ClariSure GP</b><br>Microglassfiber Cartridge Filters |
| Pre-filtration of difficult to filter SVP like Oxytetracycline  | <ul style="list-style-type: none"> <li>- High retention efficiency</li> <li>- High throughput</li> </ul>                              | -  | -   | <b>ClariSure GK</b><br>PES Membrane Cartridge Filters with Microglassfiber upstream | <b>ClariSure GP</b><br>Microglassfiber Cartridge Filters |
| Polishing of turbid solutions                                   | <ul style="list-style-type: none"> <li>- Very high retention efficiency for colloidal particles</li> <li>- High throughput</li> </ul> | <b>ClariPro GK</b><br>PES Membrane Capsule Filters with Microglassfiber upstream | <b>ClariCap GS</b><br>Microglassfiber Capsule Filters | <b>ClariSure GK</b><br>PES Membrane Cartridge Filters with Microglassfiber upstream | <b>ClariSure GS</b><br>Microglassfiber Cartridge Filters |

# ClariPro GK- 1", 2", 5", 8"

**Complies with  
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

**ClariPro GK** 0.5µm hydrophilic PES membrane capsule filters are ready to use, disposable filtration devices. These filters are specially designed filters incorporating a microglassfiber upstream layer and a downstream PES membrane layer and are used as pre- filters in biopharmaceuticals process development as well as manufacturing processes for difficult to filter solutions.

**Radiation Sterilizable:** ClariPro GK-

**Autoclavable:** ClariPro GK


## Specifications

| Construction                        |   |                     |                      |                      |
|-------------------------------------|---|---------------------|----------------------|----------------------|
| Pore Size                           | 0.5µm   |                     |                      |                      |
| Membrane                            | Hydrophilic PES   |                     |                      |                      |
| Pre-filter                          | Microglassfiber   |                     |                      |                      |
| Support Layers                      | Polyester   |                     |                      |                      |
| Body and Core                       | Polypropylene   |                     |                      |                      |
| Size                                |   |                     |                      |                      |
| Size                                | 1"  | 2"                  | 5"                   | 8"                   |
| Effective Filtration Area (Nominal) | 150 cm <sup>2</sup>   | 500 cm <sup>2</sup> | 1000 cm <sup>2</sup> | 1500 cm <sup>2</sup> |
| Vent and Drain                      | 1/4" Hose Barb with double Silicone 'O' rings   |                     |                      |                      |
| Operational                         |   |                     |                      |                      |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )  |                     |                      |                      |
| Max. Differential Pressure          | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C  |                     |                      |                      |
| By Irradiation                      | Gamma Irradiatable up to 50 kGy   |                     |                      |                      |
| By Gas                              | Sterilizable by Ethylene Oxide  |                     |                      |                      |
| Sterilization                       | Autoclavable at 121 °C for 30 minutes,<br>By Autoclave 1 cycle after gamma irradiation. Can not be in-line steam sterilized |                     |                      |                      |



## Water Flow Rates

0.5 µm ClariPro GK Capsule Filters



## Ordering Information

| Type            |      | Size |      | Pore Size |      | Inlet/Outlet             |      | Radiation Sterilizable |      | X | Sterility     |   | Pack Size |      |
|-----------------|------|------|------|-----------|------|--------------------------|------|------------------------|------|---|---------------|---|-----------|------|
|                 | Code |      | Code |           | Code |                          | Code |                        | Code |   |               |   | Code      | Code |
| ClariPro GK     | DGKX | 1"   | 51   | 0.5 µm    | 04   | ¼" SHB                   | A    | Yes                    | R    |   | Non Sterile   | 1 | 1         | 01   |
|                 |      | 2"   | 52   |           |      | ¼" MNPT                  | B    |                        |      |   | EO Sterile    | 2 |           |      |
|                 |      | 5"   | 53   |           |      | ½" MNPT                  | C    |                        |      |   | Gamma Sterile | 3 |           |      |
|                 |      | 8"   | 57   |           |      | ½" Hose Barb             | D    |                        |      |   |               |   |           |      |
|                 |      |      |      |           |      | 1½" Sanitary Flange      | E    |                        |      |   |               |   |           |      |
|                 |      |      |      |           |      | ¾" Sanitary Flange       | S    |                        |      |   |               |   |           |      |
|                 |      |      |      |           |      | Quick Connector          | J    |                        |      |   |               |   |           |      |
|                 |      |      |      |           |      | Single Step ½" Hose Barb | Q    |                        |      |   |               |   |           |      |
| <b>Example:</b> |      |      |      |           |      |                          |      |                        |      |   |               |   |           |      |
|                 | DGKX |      | 57   |           | 04   |                          | DD   |                        | R    | X |               | 2 |           | 01   |

For End Connection availability and dimensions with different sizes refer Page no. 59.

# ClariPro GK- 10", 20", 30"

**Complies with  
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

**ClariPro GK** 0.5 µm hydrophilic PES membrane large capsule filters are ready to use, disposable filtration devices with a microglassfiber upstream layer and a downstream PES membrane final layer.

The microglassfiber layer offers high dirt holding capacity and efficient retention of colloidal particles to give clear downstream in case of difficult to filter turbid solutions. The downstream PES membrane offers absolute retention, and low protein binding along with high flow rates.

**Radiation Sterilizable:** ClariPro GK-  
**Autoclavable:** ClariPro GK


## Specifications

| Construction                        |   |   |                       |
|-------------------------------------|---|---|-----------------------|
| Pore Size                           | 0.5 µm  |   |                       |
| Membrane                            | Hydrophilic PES                               |   |                       |
| Pre-filter                          | Microglassfiber                               |   |                       |
| Support Layers                      | Polyester                                     |   |                       |
| Body and Core                       | Polypropylene                                 |   |                       |
| Size                                |   |   |                       |
| Size                                | 10"   | 20"   | 30"                   |
| Effective Filtration Area (Nominal) | 5000 cm <sup>2</sup>                          | 10000 cm <sup>2</sup>   | 15000 cm <sup>2</sup> |
| Vent and Drain                      | 1/4" Hose Barb with double Silicone 'O' rings |   |                       |
| Operational                         |   |   |                       |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )      |   |                       |
| Max. Differential Pressure          | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C        |   |                       |
| Sterilization                       | By Irradiation                                | Gamma Irradiatable up to 50 kGy   |                       |
|                                     | By Gas  | Sterilizable by Ethylene Oxide  |                       |
|                                     | By Autoclave                                  | Autoclavable at 121 °C for 30 minutes, 1 cycle after gamma irradiation. Can not be in-line steam sterilized |                       |



## Water Flow Rates

0.5 µm ClariPro GK Large Capsule Filters



## Ordering Information

| Type            |      | Size |      | Pore Size |      | Inlet/Outlet             |      | Radiation Sterilizable |      | X | Inline / T-line |   | Sterility     |   | Pack Size |    |
|-----------------|------|------|------|-----------|------|--------------------------|------|------------------------|------|---|-----------------|---|---------------|---|-----------|----|
|                 | Code |      | Code |           | Code |                          | Code |                        | Code |   |                 |   | Code          |   | Code      |    |
| ClariProGK      | LGKX | 10"  | 54   | 0.5 µm    | 04   | 1½" Sanitary Flange      | E    | Yes                    | R    |   | Inline          | X | Non Sterile   | 1 | 1         | 01 |
|                 |      | 20"  | 55   |           |      | Single Step ½" Hose Barb | Q    |                        |      |   | T-line          | T | EO Sterile    | 2 |           |    |
|                 |      | 30"  | 56   |           |      |                          |      |                        |      |   |                 |   | Gamma Sterile | 3 |           |    |
| <b>Example:</b> |      |      |      |           |      |                          |      |                        |      |   |                 |   |               |   |           |    |
|                 | LGKX |      | 54   |           | 04   |                          | EE   |                        | R    |   |                 | X |               | 1 |           | 01 |

For End Connection availability and dimensions with different sizes refer Page no. 59.

# ClariSure GK

Complies with  
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds  
USFDA 21 CFR 177.1520




**ClariSure GK** 0.5 µm hydrophilic PES membrane cartridge filters offer a microglassfiber upstream layer with a PES membrane final filter to combine high dirt holding capacities with efficient retention of colloidal particles to give clear downstream in case of difficult to filter turbid solutions.

## Specifications

| Construction                        |   |
|-------------------------------------|---|
| Pore Size                           | 0.5 µm  |
| Membrane                            | Hydrophilic PES   |
| Pre-filter                          | Microglassfiber   |
| Support Layers                      | Polyester   |
| Body and Core                       | Polypropylene   |
| Size                                |   |
| Size                                | 5" 10" 20" 30"  |
| Effective Filtration Area (Nominal) | 2500 cm <sup>2</sup> 5000 cm <sup>2</sup> 10000 cm <sup>2</sup> 15000 cm <sup>2</sup> |
| Operational                         |   |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )  |
| Max. Differential Pressure          | 50 psi (3.5 Kg/cm <sup>2</sup> ) @ 25 °C  |
| Sterilization                       | Autoclavable/In-line steam sterilizable at 121 °C for 30 minutes, 25 cycles           |

## Water Flow Rates

0.5 µm ClariSure GK Cartridge Filters



## Ordering Information

| Type            |      | Size |      | Pore Size |      | Adapter |      | Elastomer              |      | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|---------|------|------------------------|------|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |         | Code |                        | Code |             | Code |           | Code |
| ClariSure GK    | CGKX | 5"   | 53   | 0.5 µm    | 04   | 7P      | A0   | Silicone               | SS   | Non Sterile | 1    | 1         | 01   |
|                 |      | 10"  | 54   |           |      | BEO     | B0   | Viton                  | SV   |             |      |           |      |
|                 |      | 20"  | 55   |           |      | 28      | C0   | EPDM                   | SE   |             |      |           |      |
|                 |      | 30"  | 56   |           |      | 'O'     | D0   | FEP Encapsulated Viton | FV*  |             |      |           |      |
|                 |      |      |      |           |      | 4463    | E0   |                        |      |             |      |           |      |
|                 |      |      |      |           |      | 4463B   | H0   |                        |      |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |         |      |                        |      |             |      |           |      |
|                 | CGKX |      | 56   |           | 04   |         | A0   |                        | SS   |             | 1    |           | 01   |

\*FV is available in Adapter Code A0 (7P) only

# ClariCap GS- 1", 2", 5", 8"

**Complies with  
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

**ClariCap GS** Capsule filters employ microglassfiber filter media for efficient retention of colloidal particles to give clear downstream in case of difficult to filter turbid solutions. These specially designed filtration devices are non media migrating with a heat calendered polypropylene layer in the downstream.




## Specifications

| Construction                        |  |                     |                      |                      |
|-------------------------------------|--|---------------------|----------------------|----------------------|
| Pore Size                           | 0.7 µm, 1 µm, 1.5 µm, 2 µm   |                     |                      |                      |
| Filter Media                        | Microglassfiber  |                     |                      |                      |
| Body and Core                       | Polypropylene  |                     |                      |                      |
| Size                                |  |                     |                      |                      |
| Size                                | 1"   | 2"                  | 5"                   | 8"                   |
| Effective Filtration Area (Nominal) | 150 cm <sup>2</sup>  | 500 cm <sup>2</sup> | 1000 cm <sup>2</sup> | 2000 cm <sup>2</sup> |
| Vent and Drain                      | 1/4" Hose Barb with double Silicone 'O' rings  |                     |                      |                      |
| Operational                         |  |                     |                      |                      |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )   |                     |                      |                      |
| Max. Differential Pressure          | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C   |                     |                      |                      |
| Sterilization                       | Autoclavable at 121 °C for 30 minutes, 30 cycles. Can not be inline steam sterilized |                     |                      |                      |

## Water Flow Rates

0.7 µm ClariCap GS Capsule Filters



## Ordering Information

| Type            |      | Size |      | Pore Size |      | Inlet/Outlet             |      | X | Bell    |      | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|--------------------------|------|---|---------|------|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |                          | Code |   |         | Code |             | Code |           | Code |
| ClariCap GS     | DGSX | 1"   | 51   | 0.7 µm    | 41   | ¼" SHB                   | A    |   | Yes     | B    | Non Sterile | 1    | 1         | 01   |
|                 |      | 2"   | 52   | 1 µm      | 05   | ¼" MNPT                  | B    |   | No Bell | X    |             |      |           |      |
|                 |      | 5"   | 53   | 1.5 µm    | 14   | ½" MNPT                  | C    |   |         |      |             |      |           |      |
|                 |      | 8"   | 57   | 2 µm      | 15   | ½" Hose Barb             | D    |   |         |      |             |      |           |      |
|                 |      |      |      |           |      | 1½" Sanitary Flange      | E    |   |         |      |             |      |           |      |
|                 |      |      |      |           |      | ¾" Sanitary Flange       | S    |   |         |      |             |      |           |      |
|                 |      |      |      |           |      | Quick Connector          | J    |   |         |      |             |      |           |      |
|                 |      |      |      |           |      | Single Step ½" Hose Barb | Q    |   |         |      |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |                          |      |   |         |      |             |      |           |      |
|                 | DGSX | 57   |      | 05        |      | DD                       |      | X | X       |      | 1           |      | 01        |      |

For End Connection availability, Bell and dimensions with different sizes refer Page no. 59.

# ClariCap GS- 10", 20", 30"

**Complies with  
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

**ClariCap GS** large capsule filters are multilayered, high throughput filters, specially designed for difficult to filter solutions.

These are high efficiency pre-filters combining the unique abilities of microglassfiber filter media to retain colloidal particles and heat calendered polypropylene filter media to ensure non media migration.




## Specifications

| Construction                        |  |                      |                      |
|-------------------------------------|--|----------------------|----------------------|
| Pore Size                           | 0.7 μm, 1 μm, 1.5 μm, 2 μm   |                      |                      |
| Filter Media                        | Microglassfiber  |                      |                      |
| Body and Core                       | Polypropylene  |                      |                      |
| Size                                |  |                      |                      |
| Size                                | 10"  | 20"                  | 30"                  |
| Effective Filtration Area (Nominal) | 3000 cm <sup>2</sup>   | 6000 cm <sup>2</sup> | 9000 cm <sup>2</sup> |
| Vent and Drain                      | 1/4" Hose Barb with double Silicone 'O' rings  |                      |                      |
| Operational                         |  |                      |                      |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )   |                      |                      |
| Max. Differential Pressure          | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C   |                      |                      |
| Sterilization                       | Autoclavable at 121 °C for 30 minutes, 30 cycles. Can not be inline steam sterilized |                      |                      |

## Water Flow Rates

0.7μm ClariCap GS Large Capsule Filters



## Ordering Information

| Type            |      | Size |      | Pore Size |      | Inlet/Outlet               |      | X | Inline / T-line |      | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|----------------------------|------|---|-----------------|------|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |                            | Code |   |                 | Code |             | Code |           | Code |
| ClariCap GS     | LGSX | 10"  | 54   | 0.7 μm    | 41   | 1 1/2" Sanitary Flange     | E    |   | Inline          | X    | Non Sterile | 1    | 1         | 01   |
|                 |      | 20"  | 55   | 1 μm      | 05   | Single Step 1/2" Hose Barb | Q    |   | T-line          | T    |             |      |           |      |
|                 |      | 30"  | 56   | 1.5 μm    | 14   |                            |      |   |                 |      |             |      |           |      |
|                 |      |      |      | 2 μm      | 15   |                            |      |   |                 |      |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |                            |      |   |                 |      |             |      |           |      |
|                 | LGSX |      | 54   |           | 41   |                            | EE   |   |                 | X    |             | 1    |           | 01   |

For End Connection availability and dimensions with different sizes refer Page no. 59.

# ClariSure GS

**Complies with  
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

**ClariSure GS** cartridge filters are multilayered, high throughput filters, specially designed for difficult to filter solutions.

It is a very high efficiency pre-filter combining the unique abilities of microglassfiber filter media to retain colloidal particles and heat calendered polypropylene filter media to ensure non media migration.




## Specifications

| Construction                        |   |
|-------------------------------------|---|
| Pore Size                           | 0.7 µm  |
| Filter Media                        | Microglassfiber   |
| Support Layers                      | Polyester   |
| Body and Core                       | Polypropylene   |
| Size                                |   |
| Size                                | 5"    10"    20"    30"   |
| Effective Filtration Area (Nominal) | 1500 cm <sup>2</sup> 3000 cm <sup>2</sup> 6000 cm <sup>2</sup> 9000 cm <sup>2</sup> |
| Operational                         |   |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )  |
| Max. Differential Pressure          | 50 psi (3.5 Kg/cm <sup>2</sup> ) @ 30 °C  |
| Sterilization                       | Autoclavable/In-line steam sterilizable at 121 °C for 30 minutes, 30 cycles         |

## Water Flow Rates

0.7 µm ClariSure GS Cartridge Filters



## Ordering Information

| Type            |      | Size |      | Pore Size |      | Adapter |      | Elastomer          |      | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|---------|------|--------------------|------|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |         | Code |                    | Code |             | Code |           | Code |
| ClariSure GS    | CPGS | 5"   | 53   | 0.7 µm    | 41   | 7P      | A0   | Silicone           | SS   | Non Sterile | 1    | 1         | 01   |
|                 |      | 10"  | 54   |           |      | BEO     | B0   | Viton              | SV   |             |      |           |      |
|                 |      | 20"  | 55   |           |      | 28      | C0   | EPDM               | SE   |             |      |           |      |
|                 |      | 30"  | 56   |           |      | 'O'     | D0   | FEP                | FV*  |             |      |           |      |
|                 |      |      |      |           |      | 4463    | E0   | Encapsulated Viton |      |             |      |           |      |
|                 |      |      |      |           |      | 4463B   | H0   |                    |      |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |         |      |                    |      |             |      |           |      |
|                 | CPGS | 55   |      | 41        |      | A0      |      | SV                 |      | 1           |      | 01        |      |

\* FV is available in Adapter Code A0 (7P) only



# ClariSure GP

**Complies with  
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

**ClariSure GP** cartridge filters are specially designed multilayered precision filtration devices. These cartridge filters act as throughput enhancers, specially in case of difficult to filter solutions.

A microglassfiber upstream layer retains very fine colloidal particles and a downstream polypropylene layer checks any kind of media migration.




## Specifications

| Construction                        |   |                      |                      |                       |
|-------------------------------------|---|----------------------|----------------------|-----------------------|
| Pore Size                           | 0.5 μm, 1 μm, 2 μm  |                      |                      |                       |
| Filter Media                        | Microglassfiber   |                      |                      |                       |
| Support Layers                      | Polypropylene   |                      |                      |                       |
| Body and Core                       | Polypropylene   |                      |                      |                       |
| Size                                |   |                      |                      |                       |
| Size                                | 5"  | 10"                  | 20"                  | 30"                   |
| Effective Filtration Area (Nominal) | 2500 cm <sup>2</sup>  | 5000 cm <sup>2</sup> | 1000 cm <sup>2</sup> | 15000 cm <sup>2</sup> |
| Operational                         |   |                      |                      |                       |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )                                    |                      |                      |                       |
| Max. Differential Pressure          | 50 psi (3.5 Kg/cm <sup>2</sup> ) @ 25 °C                                    |                      |                      |                       |
| Sterilization                       | Autoclavable/In-line steam sterilizable at 121 °C for 30 minutes, 30 cycles |                      |                      |                       |

## Water Flow Rates

ClariSure GP, 10" Cartridge Filters



## Ordering Information

| Type         |      | Size |      | Pore Size |      | Adapter |      | Elastomer              |      | Sterility   |      | Pack Size |      |
|--------------|------|------|------|-----------|------|---------|------|------------------------|------|-------------|------|-----------|------|
|              | Code |      | Code |           | Code |         | Code |                        | Code |             | Code |           | Code |
| ClariSure GP | CPGP | 5"   | 53   | 0.5 μm    | 04   | 7P      | A0   | Silicone               | SS   | Non Sterile | 1    | 1         | 01   |
|              |      | 10"  | 54   | 1 μm      | 05   | BEO     | B0   | Viton                  | SV   |             |      |           |      |
|              |      | 20"  | 55   | 2 μm      | 15   | 28      | C0   | EPDM                   | SE   |             |      |           |      |
|              |      | 30"  | 56   |           |      | 'O'     | D0   | FEP Encapsulated Viton | FV*  |             |      |           |      |
|              |      |      |      |           |      | 4463    | E0   |                        |      |             |      |           |      |
|              |      |      |      |           |      | 4463B   | H0   |                        |      |             |      |           |      |

**Example:**

|             |           |           |           |           |          |           |
|-------------|-----------|-----------|-----------|-----------|----------|-----------|
| <b>CPGP</b> | <b>55</b> | <b>04</b> | <b>A0</b> | <b>SS</b> | <b>1</b> | <b>01</b> |
|-------------|-----------|-----------|-----------|-----------|----------|-----------|

\* FV is available in Adapter Code A0 (7P) only

# ClariCap PP- 1", 2", 5", 8"

ClariCap PP capsule filters employ high retention efficiency polypropylene filter media with wide chemical compatibility for pre-filtration and polishing applications.

**Complies with  
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**




## Specifications

| Construction                        |  |                     |                      |                      |
|-------------------------------------|--|---------------------|----------------------|----------------------|
| Pore Size                           | 0.5 µm, 1 µm, 2.5 µm, 5 µm, 10 µm, 20 µm   |                     |                      |                      |
| Filter Media                        | Polypropylene  |                     |                      |                      |
| Support Layers                      | Polypropylene  |                     |                      |                      |
| Body and Core                       | Polypropylene  |                     |                      |                      |
| Size                                |  |                     |                      |                      |
| Size                                | 1"   | 2"                  | 5"                   | 8"                   |
| Effective Filtration Area (Nominal) | 250 cm <sup>2</sup>  | 500 cm <sup>2</sup> | 1000 cm <sup>2</sup> | 2000 cm <sup>2</sup> |
| Vent and Drain                      | 1/4" Hose Barb with double Silicone 'O' rings  |                     |                      |                      |
| Operational                         |  |                     |                      |                      |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )   |                     |                      |                      |
| Max. Differential Pressure          | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C   |                     |                      |                      |
| Sterilization                       | Autoclavable at 121 °C for 30 minutes, 30 cycles. Can not be inline steam sterilized |                     |                      |                      |

## Water Flow Rates

ClariCap PP, 5" Capsule Filters



## Ordering Information

| Type            |      | Size |      | Pore Size |      | Inlet/Outlet             |      | X  | Bell    |      | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|--------------------------|------|----|---------|------|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |                          | Code |    |         | Code |             | Code |           | Code |
| ClariCap PP     | DOLX | 1"   | 51   | 0.5 µm    | 04   | ¼" SHB                   | A    |    | Yes     | B    | Non Sterile | 1    | 1         | 01   |
|                 |      | 2"   | 52   | 1 µm      | 05   | ¼" MNPT                  | B    |    | No Bell | X    |             |      |           |      |
|                 |      | 5"   | 53   | 2.5 µm    | 06   | ½" MNPT                  | C    |    |         |      |             |      |           |      |
|                 |      | 8"   | 57   | 5 µm      | 07   | ½" Hose Barb             | D    |    |         |      |             |      |           |      |
|                 |      |      |      | 10 µm     | 08   | 1½" Sanitary Flange      | E    |    |         |      |             |      |           |      |
|                 |      |      |      | 20 µm     | 11   | ¾" Sanitary Flange       | S    |    |         |      |             |      |           |      |
|                 |      |      |      |           |      | Quick Connector          | J    |    |         |      |             |      |           |      |
|                 |      |      |      |           |      | Single Step ½" Hose Barb | Q    |    |         |      |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |                          |      |    |         |      |             |      |           |      |
|                 | DOLX | 53   | 05   | CC        | X    | X                        | 1    | 01 |         |      |             |      |           |      |

For End Connection availability, Bell and dimensions with different sizes refer Page no. 59.

# ClariCap PP- 10", 20", 30"

**Complies with  
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

**ClariCap PP** large capsule filters employ high retention efficiency Polypropylene filter media for wide chemical compatibility, efficient pre-filtration, with clarification and polishing applications.




## Specifications

| Construction                        |  |                       |                       |
|-------------------------------------|--|-----------------------|-----------------------|
| Pore Size                           | 0.5 µm, 1 µm, 2.5 µm, 5 µm, 10 µm, 20 µm   |                       |                       |
| Filter Media                        | Polypropylene  |                       |                       |
| Support Layers                      | Polypropylene  |                       |                       |
| Body and Core                       | Polypropylene  |                       |                       |
| Size                                |  |                       |                       |
| Size                                | 10"  | 20"                   | 30"                   |
| Effective Filtration Area (Nominal) | 5000 cm <sup>2</sup>   | 10000 cm <sup>2</sup> | 15000 cm <sup>2</sup> |
| Vent and Drain                      | 1/4" Hose Barb with double Silicone 'O' rings  |                       |                       |
| Operational                         |  |                       |                       |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )   |                       |                       |
| Max. Differential Pressure          | 60 psi (4 Kg/cm <sup>2</sup> ) @ 30 °C   |                       |                       |
| Sterilization                       | Autoclavable at 121 °C for 30 minutes, 30 cycles. Can not be inline steam sterilized |                       |                       |

## Water Flow Rates

ClariCap PP 10", Large Capsule Filters



### End Connection Type:

E: 1½" Sanitary Flange

## Ordering Information

| Type            |      | Size |      | Pore Size |      | Inlet/Outlet             |      | X  | Inline / T-line |      | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|--------------------------|------|----|-----------------|------|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |                          | Code |    |                 | Code |             | Code |           | Code |
| ClariCap PP     | LOLX | 10"  | 54   | 0.5 µm    | 04   | 1½" Sanitary Flange      | E    |    | Inline          | X    | Non Sterile | 1    | 1         | 01   |
|                 |      | 20"  | 55   | 1 µm      | 05   | Single Step ½" Hose Barb | Q    |    | T-line          | T    |             |      |           |      |
|                 |      | 30"  | 56   | 2.5 µm    | 06   |                          |      |    |                 |      |             |      |           |      |
|                 |      |      |      | 5 µm      | 07   |                          |      |    |                 |      |             |      |           |      |
|                 |      |      |      | 10 µm     | 08   |                          |      |    |                 |      |             |      |           |      |
|                 |      |      |      | 20 µm     | 11   |                          |      |    |                 |      |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |                          |      |    |                 |      |             |      |           |      |
|                 | LOLX | 54   | 06   | QQ        | X    | T                        | 1    | 01 |                 |      |             |      |           |      |

For End Connection availability and dimensions with different sizes refer Page no. 59.

# ClariSure PA

**Complies with  
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

ClariSure PA cartridge filters are 100% polypropylene, very high retention efficiency pleated cartridge filters offering large filtration area.

These filters have heat stable construction and are used as pre-filters to sterilizing membrane cartridge filters.




## Specifications

| Construction                        |  |                      |                       |                       |
|-------------------------------------|--|----------------------|-----------------------|-----------------------|
| Pore Size                           | 0.5 µm, 1 µm, 1.2 µm, 2.5 µm, 5 µm, 8 µm, 10 µm                              |                      |                       |                       |
| Filter Media                        | Polypropylene  |                      |                       |                       |
| Support Layers                      | Polypropylene  |                      |                       |                       |
| Body and Core                       | Polypropylene  |                      |                       |                       |
| Size                                |  |                      |                       |                       |
| Size                                | 5"   | 10"                  | 20"                   | 30"                   |
| Effective Filtration Area (Nominal) | 2500 cm <sup>2</sup>   | 5000 cm <sup>2</sup> | 10000 cm <sup>2</sup> | 15000 cm <sup>2</sup> |
| Operational                         |  |                      |                       |                       |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )                                     |                      |                       |                       |
| Max. Differential Pressure          | 50 psi (3.5 Kg/cm <sup>2</sup> ) @ 25 °C                                     |                      |                       |                       |
| Sterilization                       | Autoclavable/In-line steam sterilizable at 121 °C for 30 minutes, 100 cycles |                      |                       |                       |

## Water Flow Rates

ClariSure PA, 10" Cartridge Filters



## Ordering Information

| Type            |      | Size |      | Pore Size |      | Adapter |      | Elastomer              |      | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|---------|------|------------------------|------|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |         | Code |                        | Code |             | Code |           | Code |
| ClariSure PA    | CPPA | 5"   | 53   | 0.5 µm    | 04   | 7P      | A0   | Silicone               | SS   | Non Sterile | 1    | 1         | 01   |
|                 |      | 10"  | 54   | 1 µm      | 05   | BEO     | B0   | Viton                  | SV   |             |      |           |      |
|                 |      | 20"  | 55   | 1.2 µm    | 10   | 28      | C0   | EPDM                   | SE   |             |      |           |      |
|                 |      | 30"  | 56   | 2.5 µm    | 06   | 'O'     | D0   | FEP Encapsulated Viton | FV*  |             |      |           |      |
|                 |      |      |      | 5 µm      | 07   | 4463    | E0   |                        |      |             |      |           |      |
|                 |      |      |      | 8 µm      | 17   | 4463B   | H0   |                        |      |             |      |           |      |
|                 |      |      |      | 10 µm     | 08   |         |      |                        |      |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |         |      |                        |      |             |      |           |      |
|                 | CPPA | 56   | 05   | A0        | SS   | 1       | 01   |                        |      |             |      |           |      |

\*FV is available in Adapter Code A0 (7P) only

# Microglassfiber Disc Filters

## Microglassfiber Disc Filters Type - GF2

GF2 filters are high dirt holding microglassfiber disc filters specially designed for pre-filtration of solutions with high dirt load.

### Special Features

- High flow rates
- High dirt holding capacity

## Fine Microglassfiber Disc Filters - Type GFS

GFS filters are high retention efficiency fine microglassfiber disc filters.

### Special Features

- High retention efficiency

## Positively Charged Microglassfiber Disc Filters- Type GFSZ

GFSZ filters are positively charged high retention efficiency microglassfiber disc filters.

### Special Features

- Very high retention efficiency for negatively charged particles

## Ordering Information

| Type |      | Size   |      | Pore Size |      | XX | XX | Sterility   |      | Pack Size |      |
|------|------|--------|------|-----------|------|----|----|-------------|------|-----------|------|
|      | Code |        | Code |           | Code |    |    |             | Code |           | Code |
| GF2  | GF2X | 127 mm | 15   | 1.5 μm    | 14   |    |    | Non Sterile | 1    | 50        | 03   |
| GFS  | GFSX | 142 mm | 16   | 2 μm      | 15   |    |    |             |      |           |      |
| GFSZ | GFSZ | 257 mm | 17   |           |      |    |    |             |      |           |      |
|      |      | 279 mm | 18   |           |      |    |    |             |      |           |      |
|      |      | 293 mm | 19   |           |      |    |    |             |      |           |      |

**Example:**

|             |           |           |           |           |          |           |   |   |   |   |   |
|-------------|-----------|-----------|-----------|-----------|----------|-----------|---|---|---|---|---|
| ↓           | ↓         | ↓         | ↓         | ↓         | ↓        | ↓         | ↓ | ↓ | ↓ | ↓ | ↓ |
| <b>GFSX</b> | <b>18</b> | <b>15</b> | <b>XX</b> | <b>XX</b> | <b>1</b> | <b>03</b> |   |   |   |   |   |

# Filters for Polishing and Clarification

**mdi** offers a range of cartridge filters for polishing and clarification applications requiring absolute removal of particulate contaminants.

These are non media migrating, biologically and chemically inert, large area filter cartridges, offering very high (99.999%) to moderate high (99.9%) retention efficiencies to suit different applications.

## Types Available:

- >> **ClariSure PP**- Pleated polypropylene cartridge filters
- >> **ClariSure PL**- Pleated profile cartridge filters
- >> **ClariSure DP**- Pleated depth cartridge filters

## Applications

- Filtration of organic solvents in non sterile API
- Polishing filtration of Non Sterile API
- Filtration of precipitating agents
- Filtration of wash solvents for final non sterile API
- Filtration of air to dryers and micronizers
- Final wash water for Ampoule/Vial washing and Bung washing
- Bottle washing in oral formulations
- Filtration of feed water for RO plants

## Quality Assurance

These filter devices are manufactured in class 10,000 clean rooms under ISO 9001 : 2008 certified quality management systems and are validated to meet compendia and regulatory requirements.

| Assurance                          |  |
|------------------------------------|--|
| <b>Toxicity</b>                    | Passes Bioreactivity test, In Vivo, as per USP <88> for Class VI plastics                              |
| <b>Non Fiber Releasing</b>         | Passes test as per USP and comply with USFDA 21 CFR Part 210.3 (b)(6) for fiber release                |
| <b>Extractables with WFI</b>       | Passes test as per USP   |
| <b>Oxidizable Substances</b>       | Within limits as specified in USP  |
| <b>Particle Shedding</b>           | Passes USP test for particulates in injectables  |
| <b>Indirect Food Additive</b>      | All Polypropylene components meet the FDA Indirect Food Additive requirements cited in 21 CFR 177.1520 |
| <b>Good Manufacturing Practice</b> | These products are manufactured in a facility which adheres to Good Manufacturing Practices            |

## Filter Selection Chart

| Application Area   | Key Application Requirements  | Cartridge Filters  |  |  |
|--|---|--|--|--|
| Filtration of organic solvents in non sterile API          | <ul style="list-style-type: none"> <li>- High retention efficiency</li> <li>- Wide chemical compatibility</li> </ul>      | <b>ClariSure PP</b><br>Pleated Polypropylene Cartridge Filters | <b>ClariSure PL</b><br>Pleated Profile Cartridge Filters | <b>ClariSure DP</b><br>Pleated Depth Cartridge Filters |
| Polishing filtration of non sterile API                    | <ul style="list-style-type: none"> <li>- Very high retention efficiency</li> <li>- Wide chemical compatibility</li> </ul> | <b>ClariSure PP</b><br>Pleated Polypropylene Cartridge Filters | <b>ClariSure PL</b><br>Pleated Profile Cartridge Filters | -  |
| Filtration of precipitating agents                         | <ul style="list-style-type: none"> <li>- Very high retention efficiency</li> <li>- Wide chemical compatibility</li> </ul> | <b>ClariSure PP</b><br>Pleated Polypropylene Cartridge Filters | <b>ClariSure PL</b><br>Pleated Profile Cartridge Filters | -  |
| Filtration of wash solvents for final non sterile API      | <ul style="list-style-type: none"> <li>- Very high retention efficiency</li> <li>- Wide chemical compatibility</li> </ul> | <b>ClariSure PP</b><br>Pleated Polypropylene Cartridge Filters | <b>ClariSure PL</b><br>Pleated Profile Cartridge Filters | -  |
| Filtration of air to dryers and micronizers                | <ul style="list-style-type: none"> <li>- Very high retention efficiency</li> </ul>  | <b>ClariSure PP</b><br>Pleated Polypropylene Cartridge Filters | <b>ClariSure PL</b><br>Pleated Profile Cartridge Filters | -  |
| Final wash water for ampoule/vial washing and bung washing | <ul style="list-style-type: none"> <li>- Very high retention efficiency</li> <li>- Wide chemical compatibility</li> </ul> | <b>ClariSure PP</b><br>Pleated Polypropylene Cartridge Filters | <b>ClariSure PL</b><br>Pleated Profile Cartridge Filters | -  |
| Bottle washing in oral formulations                        | <ul style="list-style-type: none"> <li>- Very high retention efficiency</li> <li>- Wide chemical compatibility</li> </ul> | <b>ClariSure PP</b><br>Pleated Polypropylene Cartridge Filters | <b>ClariSure PL</b><br>Pleated Profile Cartridge Filters | -  |
| Filtration of feed water for RO plants                     | <ul style="list-style-type: none"> <li>- High dirt holding capacity</li> </ul>  | -  | -  | <b>ClariSure DP</b><br>Pleated Depth Cartridge Filters |

# ClariSure PP

**Complies with  
USFDA 21 CFR 210.3 (b)(6)**

**Meets and Exceeds  
USFDA 21 CFR 177.1520**

**ClariSure PP** cartridge filters are 100% polypropylene construction precision filtration devices. These are very high retention efficiency (**99.999%**) pleated cartridge filters offering large filtration area.

These filters are used as terminal filters in applications requiring highly efficient particulate removal.




## Specifications

| Construction                        |   |
|-------------------------------------|---|
| Pore Size                           | 0.5 µm, 1 µm, 1.2 µm, 2 µm, 2.5 µm, 5 µm<br>10 µm, 20 µm, 30 µm                       |
| Filter Media                        | Polypropylene   |
| Support Layers                      | Polypropylene   |
| Body and Core                       | Polypropylene   |
| Size                                |   |
| Size                                | 5"      10"      20"      30"   |
| Effective Filtration Area (Nominal) | 2500 cm <sup>2</sup> 5000 cm <sup>2</sup> 10000 cm <sup>2</sup> 15000 cm <sup>2</sup> |
| Operational                         |   |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )  |
| Max. Differential Pressure          | 50 psi (3.5 Kg/cm <sup>2</sup> ) @ 25 °C  |

## Water Flow Rates

ClariSure PP, 10" Cartridge Filters



## Ordering Information

| Type            |      | Size |      | Pore Size |      | Adapter |      | Elastomer              |      | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|---------|------|------------------------|------|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |         | Code |                        | Code |             | Code |           | Code |
| ClariSure PP    | CPPP | 5"   | 53   | 0.5 µm    | 04   | 7P      | A0   | Silicone               | SS   | Non Sterile | 1    | 1         | 01   |
|                 |      | 10"  | 54   | 1 µm      | 05   | BEO     | B0   | Viton                  | SV   |             |      |           |      |
|                 |      | 20"  | 55   | 1.2 µm    | 10   | 28      | C0   | EPDM                   | SE   |             |      |           |      |
|                 |      | 30"  | 56   | 2 µm      | 15   | 'O'     | D0   | PTFE                   | ST   |             |      |           |      |
|                 |      |      |      | 2.5 µm    | 06   | 4463    | E0   | FEP Encapsulated Viton | FV*  |             |      |           |      |
|                 |      |      |      | 5 µm      | 07   | 4463B   | H0   |                        |      |             |      |           |      |
|                 |      |      |      | 10 µm     | 08   |         |      |                        |      |             |      |           |      |
|                 |      |      |      | 20 µm     | 11   |         |      |                        |      |             |      |           |      |
|                 |      |      |      | 30 µm     | 19   |         |      |                        |      |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |         |      |                        |      |             |      |           |      |
|                 | CPPP | 55   | 04   | A0        | FV   | 1       | 01   |                        |      |             |      |           |      |

\* FV is available in Adapter Code A0 (7P) only



# ClariSure PL

ClariSure PL cartridge filters are all polypropylene high retention efficiency (99.99%) pleated cartridge filters offering large filtration area.

Complies with  
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds  
USFDA 21 CFR 177.1520




## Specifications

| Construction                        |  |                      |                       |                       |
|-------------------------------------|--|----------------------|-----------------------|-----------------------|
| Pore Size                           | 1 µm, 2.5 µm, 5 µm, 10 µm                |                      |                       |                       |
| Filter Media                        | Polypropylene                            |                      |                       |                       |
| Support Layers                      | Polypropylene                            |                      |                       |                       |
| Body and Core                       | Polypropylene                            |                      |                       |                       |
| Size                                |  |                      |                       |                       |
| Size                                | 5"                                       | 10"                  | 20"                   | 30"                   |
| Effective Filtration Area (Nominal) | 2500 cm <sup>2</sup>                     | 5000 cm <sup>2</sup> | 10000 cm <sup>2</sup> | 15000 cm <sup>2</sup> |
| Operational                         |  |                      |                       |                       |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> ) |                      |                       |                       |
| Max. Differential Pressure          | 50 psi (3.5 Kg/cm <sup>2</sup> ) @ 25 °C |                      |                       |                       |

## Water Flow Rates

ClariSure PL, 10" Cartridge Filters



## Ordering Information

| Type         |      | Size |      | Pore Size |      | Adapter |      | Elastomer          |      | Sterility   |      | Pack Size |      |
|--------------|------|------|------|-----------|------|---------|------|--------------------|------|-------------|------|-----------|------|
|              | Code |      | Code |           | Code |         | Code |                    | Code |             | Code |           | Code |
| ClariSure PL | CPPL | 5"   | 53   | 1 µm      | 05   | 7P      | A0   | Silicone           | SS   | Non Sterile | 1    | 1         | 01   |
|              |      | 10"  | 54   | 2.5 µm    | 06   | BEO     | B0   | Viton              | SV   |             |      |           |      |
|              |      | 20"  | 55   | 5 µm      | 07   | 28      | C0   | EPDM               | SE   |             |      |           |      |
|              |      | 30"  | 56   | 10 µm     | 08   | 'O'     | D0   | PTFE               | ST   |             |      |           |      |
|              |      |      |      |           |      | 4463    | E0   | FEP                | FV*  |             |      |           |      |
|              |      |      |      |           |      | 4463B   | H0   | Encapsulated Viton |      |             |      |           |      |

**Example:**

|      |    |    |    |    |   |    |
|------|----|----|----|----|---|----|
| CPPL | 54 | 07 | B0 | SV | 1 | 01 |
|------|----|----|----|----|---|----|

\*FV is available in Adapter Code A0 (7P) only

## ClariSure DP

Complies with  
USFDA 21 CFR 210.3 (b)(6)

Meets and Exceeds  
USFDA 21 CFR 177.1520

**ClariSure DP** cartridge filters are specially designed medium retention efficiency, multilayered, all polypropylene depth cartridge filters offering high dirt holding capacity as well as high flow rates due to their pleated configuration.

These filters are primarily used as pre-filters but can also be used as terminal filters for non-critical applications.



## Specifications

| Construction                        |  |                      |                      |                       |
|-------------------------------------|--|----------------------|----------------------|-----------------------|
| Pore Size                           | 1 µm, 1.5 µm, 3 µm, 5 µm, 10 µm, 20 µm<br>30 µm, 40 µm |                      |                      |                       |
| Filter Media                        | Polypropylene  |                      |                      |                       |
| Support Layers                      | Polypropylene  |                      |                      |                       |
| Body and Core                       | Polypropylene  |                      |                      |                       |
| Size                                |  |                      |                      |                       |
| Size                                | 5"   | 10"                  | 20"                  | 30"                   |
| Effective Filtration Area (Nominal) | 2000 cm <sup>2</sup>                                   | 4000 cm <sup>2</sup> | 8000 cm <sup>2</sup> | 12000 cm <sup>2</sup> |
| Operational                         |  |                      |                      |                       |
| Max. Operating Temperature          | 80 °C @ ≤ 30 psi (2 Kg/cm <sup>2</sup> )               |                      |                      |                       |
| Max. Differential Pressure          | 50 psi (3.5 Kg/cm <sup>2</sup> ) @ 25 °C               |                      |                      |                       |

## Ordering Information

| Type            |      | Size |      | Pore Size |      | Adapter |      | Elastomer          |      | Sterility   |      | Pack Size |      |
|-----------------|------|------|------|-----------|------|---------|------|--------------------|------|-------------|------|-----------|------|
|                 | Code |      | Code |           | Code |         | Code |                    | Code |             | Code |           | Code |
| ClariSure DP    | CPDP | 5"   | 53   | 1 µm      | 05   | 7P      | A0   | Silicone           | SS   | Non Sterile | 1    | 1         | 01   |
|                 |      | 10"  | 54   | 1.5 µm    | 14   | BEO     | B0   | Viton              | SV   |             |      |           |      |
|                 |      | 20"  | 55   | 3 µm      | 16   | 28      | C0   | EPDM               | SE   |             |      |           |      |
|                 |      | 30"  | 56   | 5 µm      | 07   | 'O'     | D0   | PTFE               | ST   |             |      |           |      |
|                 |      |      |      | 10 µm     | 08   | 4463    | E0   | FEP                | FV*  |             |      |           |      |
|                 |      |      |      | 20 µm     | 11   | 4463B   | H0   | Encapsulated Viton |      |             |      |           |      |
|                 |      |      |      | 30 µm     | 19   |         |      |                    |      |             |      |           |      |
|                 |      |      |      | 40 µm     | 20   |         |      |                    |      |             |      |           |      |
| <b>Example:</b> |      |      |      |           |      |         |      |                    |      |             |      |           |      |
|                 | CPDP |      | 56   |           | 14   |         | A0   |                    | SS   |             | 1    |           | 01   |

\*FV is available in Adapter Code A0 (7P) only

# Chemical Compatibility

Table below shows the chemical compatibility of various process filtration products with some commonly used solvents. All products were exposed to specified chemicals for 72 hours at 25 °C. Chemical compatibility data on specific reagents is available on request.

| Reagents                      | AseptiSure<br>Cartridge Filters |    |    |       | ClariSure<br>Cartridge Filters |    |    |    |    |    | 'O' Rings/Gasket Seals |       |    |                              |
|-------------------------------|---------------------------------|----|----|-------|--------------------------------|----|----|----|----|----|------------------------|-------|----|------------------------------|
|                               | HS/KS                           | KR | NS | TH/TF | GS                             | GP | PA | PP | PL | DP | Silicone               | Viton | EP | FEP<br>Encapsulated<br>Viton |
| <b>Solvents</b>               |                                 |    |    |       |                                |    |    |    |    |    |                        |       |    |                              |
| Acetone                       | N                               | N  | G  | G     | G                              | G  | G  | G  | G  | G  | N                      | N     | G  | G                            |
| Acetonitrile                  | G                               | G  | G  | G     | G                              | G  | G  | G  | G  | G  | G                      | N     | G  | G                            |
| Benzene                       | G                               | G  | G  | G     | F                              | F  | F  | F  | F  | F  | N                      | G     | N  | G                            |
| Benzyl Alcohol                | N                               | N  | G  | G     | G                              | G  | G  | G  | G  | G  | G                      | G     | G  | G                            |
| Benzyl Alcohol 4%             | G                               | G  | G  | G     | G                              | G  | G  | G  | G  | G  | G                      | G     | G  | G                            |
| Diethyl Ether                 | G                               | G  | G  | G     | N                              | N  | N  | N  | N  | N  | N                      | N     | N  | G                            |
| Dimethylformamide             | N                               | N  | G  | G     | G                              | G  | G  | G  | G  | G  | G                      | N     | N  | G                            |
| Ethyl Acetate                 | G                               | G  | G  | G     | G                              | G  | G  | G  | G  | G  | N                      | N     | G  | G                            |
| Ethylene Glycol               | G                               | G  | G  | G     | G                              | G  | G  | G  | G  | G  | G                      | G     | G  | G                            |
| Hexane                        | G                               | G  | G  | G     | F                              | F  | F  | F  | F  | F  | N                      | G     | N  | G                            |
| Iso Propyl Alcohol            | G                               | G  | G  | G     | G                              | G  | G  | G  | G  | G  | G                      | G     | G  | G                            |
| Methanol                      | G                               | G  | G  | G     | G                              | G  | G  | G  | G  | G  | G                      | N     | G  | G                            |
| Methylene Chloride            | N                               | N  | N  | N     | N                              | N  | N  | N  | N  | N  | N                      | N     | N  | G                            |
| n-Butanol                     | G                               | G  | G  | G     | G                              | G  | G  | G  | G  | G  | N                      | G     | G  | G                            |
| Peanut oil                    | G                               | G  | G  | G     | G                              | G  | G  | G  | G  | G  | G                      | G     | G  | G                            |
| Tetrahydrofuran/Water (50:50) | N                               | N  | N  | G     | N                              | N  | N  | N  | N  | N  | N                      | N     | N  | G                            |
| Toluene                       | G                               | G  | G  | G     | G                              | G  | G  | G  | G  | G  | N                      | G     | N  | G                            |
| Trichloroethylene             | N                               | N  | N  | N     | N                              | N  | N  | N  | N  | N  | N                      | G     | N  | G                            |
| <b>Acids</b>                  |                                 |    |    |       |                                |    |    |    |    |    |                        |       |    |                              |
| Hydrochloric Acid 37%         | G                               | G  | N  | G     | G                              | G  | G  | G  | G  | G  | N                      | G     | N  | G                            |
| Hydrofluoric Acid 10%         | G                               | G  | F  | G     | N                              | N  | G  | G  | G  | G  | N                      | G     | F  | G                            |
| Nitric Acid 67%               | N                               | N  | N  | G     | G                              | G  | G  | G  | G  | G  | N                      | G     | N  | G                            |
| Nitric Acid 7%                | G                               | G  | N  | G     | G                              | G  | G  | G  | G  | G  | G                      | G     | G  | G                            |
| Sulphuric Acid 10%            | G                               | G  | F  | G     | G                              | G  | G  | G  | G  | G  | G                      | G     | G  | G                            |
| <b>Bases</b>                  |                                 |    |    |       |                                |    |    |    |    |    |                        |       |    |                              |
| Ammonium Hydroxide 25%        | N                               | G  | N  | G     | G                              | G  | G  | G  | G  | G  | G                      | G     | G  | G                            |
| Sodium Hydroxide 32%          | N                               | G  | N  | G     | G                              | G  | G  | G  | G  | G  | G                      | G     | G  | G                            |
| Potassium Hydroxide 32%       | N                               | G  | N  | G     | G                              | G  | G  | G  | G  | G  | G                      | G     | G  | G                            |

G: Good    F: Fair    N: Not recommended

# Chemical Compatibility

| Reagents                      | AseptiCap Capsule Filters |    |       | AseptiVent Capsule Filters | ClariCap Capsule Filters |    | Membrane Disc Filters | Microglassfiber Pre-filter Discs |
|-------------------------------|---------------------------|----|-------|----------------------------|--------------------------|----|-----------------------|----------------------------------|
|                               | KL/KS                     | KO | NL/NS | TF                         | GS                       | PP | NN                    | GF                               |
| <b>Solvents</b>               |                           |    |       |                            |                          |    |                       |                                  |
| Acetone                       | N                         | N  | G     | G                          | G                        | G  | G                     | G                                |
| Acetonitrile                  | G                         | G  | G     | G                          | G                        | G  | G                     | G                                |
| Benzene                       | G                         | G  | G     | G                          | F                        | F  | G                     | G                                |
| Benzyl Alcohol                | N                         | N  | G     | G                          | G                        | G  | G                     | G                                |
| Benzyl Alcohol 4%             | G                         | G  | G     | G                          | G                        | G  | G                     | G                                |
| Diethyl Ether                 | G                         | G  | G     | G                          | N                        | N  | G                     | G                                |
| Dimethylformamide             | N                         | N  | G     | G                          | G                        | G  | G                     | G                                |
| Ethyl Acetate                 | G                         | G  | G     | G                          | G                        | G  | G                     | G                                |
| Ethylene Glycol               | G                         | G  | G     | G                          | G                        | G  | G                     | G                                |
| Hexane                        | G                         | G  | G     | G                          | F                        | F  | G                     | G                                |
| Iso Propyl Alcohol            | G                         | G  | G     | G                          | G                        | G  | G                     | G                                |
| Methanol                      | G                         | G  | G     | G                          | G                        | G  | G                     | G                                |
| Methylene Chloride            | N                         | N  | N     | N                          | N                        | N  | G                     | G                                |
| n-Butanol                     | G                         | G  | G     | G                          | G                        | G  | G                     | G                                |
| Peanut oil                    | G                         | G  | G     | G                          | G                        | G  | G                     | G                                |
| Tetrahydrofuran/Water (50:50) | N                         | N  | N     | G                          | N                        | N  | G                     | G                                |
| Toluene                       | G                         | G  | G     | G                          | G                        | G  | G                     | G                                |
| Trichloroethylene             | N                         | N  | N     | N                          | N                        | N  | G                     | G                                |
| <b>Acids</b>                  |                           |    |       |                            |                          |    |                       |                                  |
| Hydrochloric Acid 37%         | G                         | G  | N     | G                          | G                        | G  | N                     | G                                |
| Hydrofluoric Acid 10%         | G                         | G  | G     | G                          | N                        | G  | G                     | N                                |
| Nitric Acid 67%               | N                         | N  | N     | G                          | G                        | G  | N                     | G                                |
| Nitric Acid 7%                | G                         | G  | G     | G                          | G                        | G  | G                     | G                                |
| Sulphuric Acid 10%            | G                         | G  | G     | G                          | G                        | G  | G                     | G                                |
| <b>Bases</b>                  |                           |    |       |                            |                          |    |                       |                                  |
| Ammonium Hydroxide 25%        | N                         | G  | N     | G                          | G                        | G  | G                     | G                                |
| Sodium Hydroxide 32%          | N                         | G  | N     | G                          | G                        | G  | G                     | F                                |
| Potassium Hydroxide 32%       | N                         | G  | N     | G                          | G                        | G  | N                     | G                                |

G: Good    F: Fair    N: Not recommended

## End Connections Availability Chart for Capsule Filters

| End Connections          | Small Capsule Filters |      |    |    |    |    | Large Capsule Filters |     |     |        |     |     |
|--------------------------|-----------------------|------|----|----|----|----|-----------------------|-----|-----|--------|-----|-----|
|                          |                       |      |    |    |    |    | Inline                |     |     | T-line |     |     |
|                          | 37mm                  | 50mm | 1" | 2" | 5" | 8" | 10"                   | 20" | 30" | 10"    | 20" | 30" |
| ½" Hose Barb             | X                     | X    | √  | √  | √  | √  | X                     | X   | X   | X      | X   | X   |
| Single Step ½" Hose Barb | X                     | X    | X  | √  | √  | √  | √                     | √   | √   | X      | X   | X   |
| ¼" Stepped Hose Barb     | √                     | √    | √  | √  | √  | √  | X                     | X   | X   | X      | X   | X   |
| 1½" Sanitary Flange      | X                     | X    | √  | √  | √  | √  | √                     | √   | √   | √      | √   | √   |
| ¾" Sanitary Flange       | X                     | √    | X  | √  | √  | √  | X                     | X   | X   | X      | X   | X   |
| ½" MNPT                  | X                     | X    | X  | √  | √  | √  | X                     | X   | X   | X      | X   | X   |
| ¼" MNPT                  | X                     | X    | √  | √  | √  | √  | X                     | X   | X   | X      | X   | X   |
| 1/8" MNPT                | X                     | √    | X  | X  | X  | X  | X                     | X   | X   | X      | X   | X   |
| Quick Connector          | X                     | X    | √  | √  | √  | √  | X                     | X   | X   | X      | X   | X   |

**Bell at Outlet is available with only 1" Capsule Filter with ¼" SHB End Connection**

## Dimensions: Capsule Filters

| End Connections                                   | Capsule Filter Size |        |         |       |        |        |        |        |
|---|---------------------|--------|---------|-------|--------|--------|--------|--------|
|   | 25mm                | 37mm   | 50mm    | 1"    | 2"     | 5"     | 8"     | 10"    |
| ¼" SHB I/O  | -                   | 64 mm  | 64.7 mm | 94 mm | 121 mm | 173 mm | 231 mm | -      |
| ¾" Sanitary Flange Inlet I/O                      | -                   | -      | 51.4 mm | 91 mm | 113 mm | 164 mm | 215 mm | -      |
| Female Luer Lock Inlet/<br>Male Luer Slip Out let | 23 mm               | -      | -       | -     | -      | -      | -      | -      |
| 1½" Sanitary Flange I/O                           | -                   | -      | -       | 93 mm | 111 mm | 161 mm | 213 mm | 332mm  |
| ½" Hose Barb I/O                                  | -                   | -      | -       | 91 mm | 113 mm | 164 mm | 215 mm | -      |
| ½" Single Step Hose Barb I/O                      | -                   | -      | -       | 91 mm | 113 mm | 164 mm | 215 mm | 337 mm |
| 1½" Sanitary Flange Inlet<br>½" Hose Barb Outlet  | -                   | -      | -       | 96 mm | 118 mm | 175 mm | 228 mm | 332 mm |
| Operational Radius                                | 15mm                | 22.5mm | 28mm    | 30 mm | 65 mm  | 65 mm  | 65 mm  | 78 mm  |

## Ordering Information

### **Shipment details for customers outside India**

Through Federal Express, UPS, or DHL courier (specify complete street address).

By air freight for large quantities (specify airport of discharge).

Goods usually reach destination within 5-10 days from date of shipment.

Membrane products are light weight and air freight charges usually vary between 3% to 10% of the value. Any duties/taxes in the country of destination are the responsibility of the consignee.

### **Shipment details for customers in India**

The consignments can be sent through courier. Courier charges will be borne by the customer. Please specify the preferred courier and provide any form and instructions for octroi etc. that may be required for shipment.

### **How to order**

Orders may be placed by phone/fax/email/mail directly to Sales.

### **Advanced Microdevices Pvt. Ltd.**

20-21, Industrial Area, Ambala Cantt - 133 006, INDIA

**Tel:** +91-171-2699290, 2699471

**Fax:** +91-171-2699221, 2699008

**Email:** [orders@mdimembrane.com](mailto:orders@mdimembrane.com)

[support@mdimembrane.com](mailto:support@mdimembrane.com)

## Quality Policy

Quality is built into **mdi** products and services by not only adhering to well designed quality systems to consistently produce high quality, internationally acceptable products but also by striving to incorporate superior performance parameters into all our products and services and provide our customers with a unique performance advantage in their application. Our quality policy provides a glimpse of our commitment:

"**mdi** strives to provide to its customers products and services of highest standards possible, consistently superior, and more satisfying than competing products and complying with quality management systems."



## Stride Towards Excellence

At **mdi**, our mission is to constantly strive to achieve excellence in all our endeavors by establishing systems to create excellent products and services to fulfil the needs of our customers. To achieve this we

- Frequently compare our products with competing brands
- Simulate tests for functional use
- Develop easy-to-use innovative products

We are constantly working on improvements and welcome suggestions from our customers.

## Guarantee

All **mdi** products are guaranteed and are backed by our

- Technical expertise and experience of over 30 years
- 'Special **mdi** process' for consistency and repeatability
- Strict quality control and quality assurance regimen
- Certificate of Analysis accompanying all shipments

We have an unconditional replacement policy in case of any defects.





## WORLD WIDE EXPORTS

### **MUMBAI**

Tel/ Fax : 022-28654435  
022-40214436  
Mobile : 09323801794

### **Delhi**

Mobile : 08295618833

### **AHMEDABAD**

Mobile : 09328257987

### **BANGALORE**

Mobile : 09972587761

### **CHENNAI**

Mobile : 09972587761

### **HYDERABAD**

Mobile : 09391935423

## **Advanced Microdevices Pvt. Ltd.**

20-21, Industrial Area, Ambala Cantt-133 006, INDIA  
Tel : +91-171-2699290, 2699471 Fax : +91-171-2699221, 2699008  
E-mail : [info@mdimembrane.com](mailto:info@mdimembrane.com)  
Internet : [www.mdimembrane.com](http://www.mdimembrane.com)