COLUMBIA AGAR (Sheep blood 5%)

Medium for fastidious microorganisms isolation.

TYPICAL FORMULA (g/l)

Pancreatic digest of Casein	12.0
Starch	1.0
Peptic digest of Animal tissue	5.0
Sodium Chloride	5.0
Yeast Extract	3.0
Beef extract	3.0
Sheep blood defibrinated	50.0 ml
Agar	14.0
Final pH 7.3 ± 0.2	

DESCRIPTION

COLUMBIA AGAR (Sheep blood 5%) is a medium for fastidious microorganisms isolation.

PRINCIPLE

Pancreatic digest of Casein, beef extract and Peptic digest of Animal tissue provide nitrogen, carbon, sulfur and other essential growth factors. Sodium chloride maintains the osmotic balance of the medium. Sheep blood defibrinated supplies additional growth factors for fastidious microorganisms and allows to evidence the haemolytic reactions. Agar is the solidifying agent.

TECHNIQUE

Inoculate the plates by streaking the sample to examine onto the surface of the medium using a sterile loop in order to isolate single colonies. Then stab the agar several times in order to depositate beta-haemolytic streptococci below the agar surface. Incubate the plates at 36±1°C for 18-48 hours, in aerobic, anaerobic or microaerophilic atmosphere according to the established laboratory procedure.

INTERPRETATION OF RESULTS

Observe for growth and haemolytic reactions. Four different kinds of haemolysis can be distinguished:

- alfa-haemolysis: haemoglobin is reduced to metahaemoglobin in the medium that surrounds the colony with a consequent greenish decolouring of the medium;
- beta-haemolysis: it is the lysis of erythrocytes which is evident in a light zone around the colony;
- gamma-haemolysis: it does not occur any destruction of erythrocytes or any change in the medium;
- alfa ' -haemolysis: it's evident a small zone of complete haemolysis surrounded by an area of partial lysis

STORAGE

2-8°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product is not classified as hazardous by current legislation and does not contain harmful substances in concentrations of ≥1%. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL of WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

1. Ellner, P.D., C.J. Stoessel., E. Drakeford, and F. Vasi (1966). A new culture medium for medical bacteriology. Am. J.Clin. Path. 45, 502-504.
2. Isenberg, H.D. (ed.) (1992). Clinical microbiology procedures handbook, vol. 1 American Society for Microbiology, Washington, DC.

PRODUCT SPECIFICATIONS

NAME

COLUMBIA AGAR (Sheep blood 5%)

PRESENTATION

Ready plates (90 mm) containing 20+/-1 ml of medium.

STORAGE

2-8°C

PACKAGING

Code	Content	Packaging
11025	20 plates	10 plates in thermically soldered film 2 x 10 plates in cardboard boxes
11025*	100 plates	10 plates in thermically soldered film 10 piles in box

pH OF THE MEDIUM

7.3 ± 0.2

USE

COLUMBIA AGAR (Sheep blood 5%) is a medium for fastidious microorganisms isolation.

TECHNIQUE

Refer to technical sheet of the product.

APPEARANCE OF THE MEDIUM

Cherry red medium, opaque.

SHELF LIFE

90 days

QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
7 days at 25 ± 1°C, in aerobiosis
7 days at 36 ± 1°C, in aerobiosis
- Microbiological control
Inoculum for productivity: 10-100 UFC/ml
Inoculum for selectivity: 10⁴-10⁵ UFC/ml
Inoculum for specificity: ≤ 10⁴ UFC/ml
Incubation conditions: 18-24 h at 36 ± 1°C

Microorganism		Growth	Haemolysis
<i>Staphylococcus aureus</i>	ATCC 25923	good	β-haemolysis
<i>Streptococcus pneumoniae</i>	ATCC 6305	good	α-haemolysis
<i>Streptococcus pyogenes</i>	ATCC 19615	good	β-haemolysis
<i>Escherichia coli</i>	ATCC 25922	good	

TABLE OF SYMBOLS

LOT	Batch code		Temperature limitation		Manufacturer		Contains sufficient for <n> tests
REF	Catalogue number		Keep away from heat		Use by		Caution, consult accompanying documents

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