

StarLine™ Immunoaffinity Columns

Romer Labs® offers a wide range of immunoaffinity columns for regulated mycotoxins. StarLine™ immunoaffinity columns are designed for a simple and reliable clean up before detection of mycotoxins in a broad range of foodstuffs and animal feed. StarLine™ immunoaffinity columns have been validated in an extensive testing program and can be used for the purification and concentration of mycotoxins prior to analyses by various techniques, such as HPLC, LC-MS, GC-MS, ELISA or direct fluorometry.



Mycotoxin Immunoaffinity Columns Portfolio

3 ml Format	
AflaStar™ R	Aflatoxin B1, B2, G1, G2
AflaStar™ M1 R	Aflatoxin M1
OchraStar™ R	Ochratoxin A
FumoniStar™	Fumonisin B1, B2, B3
ZearaStar™ R	Zearalenone
DONStar™ R	Deoxynivalenol
T2/HT2Star™ R	T-2 Toxin, HT-2 Toxin
1 ml Format	
AflaStar™ Fit	Aflatoxin B1, B2, G1, G2
OchraStar™ Fit	Ochratoxin A

Why use StarLine™ immunoaffinity columns?

Superior Quality

Exceed international standards and performance criteria

Accurate

Low % CV and detection limit

Reliable

Use specific monoclonal antibodies

Durable

Long shelf life

Flexible

Can be used with different extraction solutions and for different contamination levels

Versatile

For a broad range of commodities and foodstuffs

Convenient


Various application notes available



StarLine™ Immunoaffinity Columns
MYCOTOXINS

Procedure of Romer Labs® Immunoaffinity Columns

3-step clean up – StarLine™ principle



Romer Labs® Mycotoxin Immunoaffinity Columns:

Item Name	Order Number	Analytes	Format	Package Size
AflaStar™ Fit	COIAC1001	Aflatoxin B1, B2, G1, G2	1 mL	25 columns
AflaStar™ R	COIAC1004	Aflatoxin B1, B2, G1, G2	3 mL	25 columns
AflaStar™ M1 R	COIAC1005	Aflatoxin M1	3 mL	25 columns
AflaStar™ Fit	COIAC1501	Aflatoxin B1, B2, G1, G2	1 mL	500 columns
AflaStar™ R	COIAC1504	Aflatoxin B1, B2, G1, G2	3 mL	500 columns
OchraStar™ Fit	COIAC2001	Ochratoxin A	1 mL	25 columns
OchraStar™ R	COIAC2004	Ochratoxin A	3 mL	25 columns
FumoniStar™	COIAC3000	Fumonisin B1, B2 and B3	3 mL	25 columns
ZearaStar™ R	COIAC4004	Zearalenone	3 mL	25 columns
DONStar™ R	COIAC5004	Deoxynivalenol	3 mL	25 columns
DONStar™ R	COIAC5504	Deoxynivalenol	3 mL	500 columns
T2/HT2Star™ R	COIAC6004	T-2 Toxin and HT-2 Toxin	3 mL	25 columns

For further information, please contact:

• office@romerlabs.com • salesasia@romerlabs.com • office-europe@romerlabs.com