

ESC Svab

ESC (Easy Surface Checking) svab ürün grubu
gıda ve farmasötik endüstrilerinde yüzeylerin
mikrobiyolojik kontrolüne yönelik yeni nesil
örnek alma çözümleri sunmaktadır.



Ürün özellikleri

- Diluent içeren tüp ve kapağa tutturulmuş svab
- Rayon svab ucu
- Geniş diluent seçeneği
- Sızdırmaz kapak
- Oda sıcaklığında muhafaza edilebilir
- Uluslararası standartlara uygundur (ISO 18593, ISO 17604)

Ref.	Diluent	Ambalaj
85601	Neutralizing Rinse Solution	100 x 10 mL
85602	Neutralizing Rinse Solution	100 x 5 mL
85603	Buffered Peptone Water	100 x 10 mL
85604	Buffered Peptone Water	100 x 5 mL
85605	D/E Neutralizing Broth	100 x 10 mL
85606	D/E Neutralizing Broth	100 x 5 mL
85607	Lethen Broth	100 x 10 mL
85608	Lethen Broth	100 x 5 mL
85609	Maximum Recovery Diluent	100 x 10 mL
85610	Maximum Recovery Diluent	100 x 5 mL

Referanslar

1. ISO 18593:2004 Microbiology of food and animal feeding stuffs - Horizontal methods for sampling techniques from surfaces using contact plates and swabs. International Standards Organisation, Geneva
2. ISO 17604:2003 Microbiology of food and animal feeding stuffs - Carcass sampling for microbiological analysis. International Standards Organisation, Geneva
3. ISO 6887-1:1999 Microbiology of food and animal feeding stuffs - Preparation of test samples, initial suspension and decimal dilution for microbiological examination. International Standards Organisation, Geneva

